


Express Mail Label No.: EV475172226US

Date of Deposit: February 9, 2005

Attorney Docket No. 24299-502



IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANTS: Duff *et al.*

APPLICATION NUMBER: 09/247,874

EXAMINER: Schnizer, Richard A.

FILING DATE: February 10, 1999

ART UNIT: 1635

FOR: THERAPEUTICS AND DIAGNOSTICS BASED ON A NOVEL IL-1B MUTATION

MAIL STOP AF

Commissioner for Patents

P.O. Box 1450

Alexandria, VA 22313-1450

DECLARATION OF DR. FRANCESCO S. DI GIOVINE UNDER 37 CFR 1.132

Sir:

I, Francesco S. di Giovine, born on June 4, 1956, a University Senior Lecturer at the University of Sheffield, do hereby declare that:

1. I am one of the inventors of the above-identified application entitled "Therapeutics and Diagnostics Based on a Novel IL-1B Mutation."
2. I have an M.D. degree in Medicine and Surgery from the University of Florence (Italy), Faculty of Medicine, 1982, and a PhD degree in Molecular Immunology from the University of Edinburgh (UK), 1988, Faculty of Medicine.
3. I am actively engaged in researching genetic predispositions to various inflammatory diseases.
4. Working under my direction, members of my laboratory in early 1996 discovered and sequenced the IL-1B allele having a "C" rather than a "G" at the position corresponding to +6912 in Figure 1, which is nucleotide 8845 of SEQ ID NO: 1. This allele is also referred to as IL-1B (+6912) allele 2.
5. A copy of the laboratory notebook ("the Notebook") kept by my assistant Carol Campbell and reviewed by me on an ongoing basis is attached hereto as Appendix A. This notebook contains experiments performed in my laboratory from August 24, 1995 to March 22, 1996. The page numbering is located in the lower left or right-hand corner of each sheet. I confirm that this is a true and complete copy of the notebook kept by Ms. Campbell during the time-period in question.

Applicants: Duff *et al.*
U.S.S.N.: 09/247,874

6. As described in the instant application (Example 1, pp. 36-37), a PCR product corresponding to the 3'UTR of the IL-1B gene was amplified from human genomic DNA and sequenced. The primers used in the amplification were disclosed in the application as SEQ ID NO: 3 (5'-GTCCCCACATTCTGATGAGCAAC-3') and SEQ ID NO: 4 (5'-TGCAGCACTCAGCAATGAGGAG-3'), which bind to regions of the IL-1B gene corresponding to positions +6720 to +6742 and +7102 to +7123, when the IL-1B gene is numbered in accordance with the numbering of Figure 1. These primers were designed on August 29, 1995 and are represented as primers F₂ and B₁ on page 3 of the Notebook. The determination that a G to C change at the +6912 location was performed by sequencing of the fragment amplified by the F₂ and B₁ PCR primers, the result of which was recorded on February 21, 1996 on page 116 of the Notebook. The oligonucleotide primer used for sequencing of the +6912 allele 2 is shown on page 117. A sequence-specific oligonucleotide primer hybridizing to nucleotides +6913 to +6947 of the IL-1B gene is described on page 118. The location of the G to C change is further evidenced on page 124, entitled "Further plans for the +8845 polymorphism," because as stated above, nucleotide 8845 of SEQ ID NO: 1 corresponds to the IL-1B (+6912) location.

7. I understand that the Examiner, while admitting that we have discovered a G to C polymorphism at position +6912, stating that neither the specification or my previously filed declaration supports a 9721 nucleotide sequence of IL-1B with a C at position +6912. (See, Office action at page 5). The Examiner states that "[s]ince polymorphisms can occur throughout a molecule, one cannot assume that there are no other polymorphisms linked to position +6912 within the 9721 bases of the IL-1B gene, and that the sequence of the rest of the 9721 nucleotides is identical to that reported in the prior art." (See, Office action at pages 5-6).

8. I believe that the C polymorphism at position +6912 is a single nucleotide polymorphism, or "SNP." A SNP is known in the art as a DNA sequence variation among individuals in which the purine or pyrimidine base (as guanine) of a single nucleotide in the genome has been replaced by another such base (as cytosine). Therefore, I believe that one of skill in the art would recognize that the sequence of the rest of the 9721 nucleotides of IL-1B is identical to the wild-type sequence, which has been reported in the art, and that upon identifying the SNP at position +6912, it was unnecessary to re-sequence the entire IL-1B gene.

9. Regarding the Examiner's statement regarding linkage analysis of the IL-1B +6912 with other IL-1B polymorphisms, I note that genetic linkage between alleles of a given gene such as IL-1B does not indicate the presence of multiple sequence variations in the IL-1B gene of an individual. Rather, I believe that linkage analysis is useful to demonstrate co-segregating polymorphisms that contribute to a given disease or disorder. As shown in the Notebook, under my direction, linkage analysis was performed on or about February 23, 1996 between the IL-1B (+6912) locus and the IL-1B *tag* locus, and the result is shown on page 129 of the Notebook. This analysis demonstrates that the IL-1B (+6912) allele 1 (termed "Allele (G)" in the notebook) is 100% associated with the IL-1B *tag* allele 2, while the IL-1B (+6912) allele 2 (termed "Allele (C)" in the notebook) is 100% associated with the IL-1B *tag* allele 1, further

Applicants: Duff *et al.*
U.S.S.N.: 09/247,874

demonstrating that my laboratory had identified a novel allele having a "C" at position +6912 of the IL-1B gene.

10. At the time of these experiments, the human IL-1B gene sequence published by Clark *et al.* (Nucleic Acids Research 14(20):7897-7914 (1986)) was regarded as the standard sequence for human IL-1B. This sequence is also deposited in GenBank under the accession number X04500. The Clark *et al.* sequence shows a "G" at position +6912, in contrast to the allele we discovered, which contains a "C" at position +6912. We named the "G" variation "allele 1" and the "C" variation "allele 2." The statement in the application indicating that the IL-1B allele 1 has a cytosine at position +6912 and that allele 2 has a guanine at that position is a typographical error. I believe that one of skill in the art would recognize the existence of this error based on the teachings of the application and, further, in view of the contents of the Notebook. Our measurements of allele frequency presented in the patent application (*e.g.*, Example 2, pp. 37-38) demonstrate that allele 1 is the more frequent allele and may therefore be considered the wild-type allele. Moreover, multiple publications and database entries have presented the Clark *et al.* nucleic acid sequence as the wild-type sequence. (See, *e.g.*, US Patent numbers 5,686,246; 6,720,141; 6,730,476; and 6,746,839; and GenBank Accession number P01584).

11. In conclusion, in February, 1996, prior to the time of filing of the present application, I was in possession of a novel IL-1B sequence with a "C" at position +6912 from the transcription start site, which is position 8845 of a nucleic acid sequence numbered in accordance with SEQ ID NO: 1. This allele is termed the IL-1B (+6912) allele 2 and is substantially less common than the IL-1B (+6912) allele 1 known in the prior art.

Applicants: Duff *et al.*
U.S.S.N.: 09/247,874

12. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

Francesco di' Giovine

FRANCESCO S. DI GIOVINE, M.D., Ph.D.

Dated: January 31st, 2005

TRA 1988391v1

Graham

2855-322

Perkin elmer. Seq

User No
13

Liob.

PCR lab.

23

Please refer to your Niceday catalogue for full details of the Niceday Manuscript Book range

Mohammed - Work 0116 254 1414 beep 4099
- Home 0116 232 2675

Franco Hotel Fax:- 03 9622 8877

10 Mrs diff - Fax @

12pm → 10am

Niceday Manuscript Books

Range of Sizes and Formats

Size	Format	No. of Sheets	Product Code
Everyday Range			
A6	Feint & indexed A-Z	96	No.237204
A5	Feint & indexed A-Z	96	No.217539
A4	Feint & indexed A-Z	96	No.215221
A6	Feint Ruled	96	No.118080
A5	Feint Ruled	96	No.144300
A4	Feint Ruled	96	No.142610
A6	Feint & indexed A-Z Wirebound	80	No.287533
A5	Feint & indexed A-Z Wirebound	80	No.178695
A4	Feint & indexed A-Z Wirebound	80	No.164569
A6	Feint Ruled Wirebound	80	No.199230
A5	Feint Ruled Wirebound	80	No.179273
A4	Feint Ruled Wirebound	80	No.233139

Environment Range (100% Recycled)

A6	Feint Ruled	96	No.174951
A5	Feint Ruled	96	No.283590
A4	Feint Ruled	96	No.118367

Economy Range

A6	Feint Ruled	96	No.135969
A5	Feint Ruled	96	No.174741
A4	Feint Ruled	96	No.240288

GRANT No. [56833]

use 100g/m
w/ 500000
w/ 500000
need 100g/m
w/ 500000
w/ 500000

(M) PNAS - 1993 90 6 2295.

(D) J. o. immun 94 153-2 712.

(M) ADVANCES IN NEUROIMMUNOLOGY 92 2 No 1 -


(D) Int J. immuno. pharm. 42 16 No 3

(M) MOL & CEL BIOG 95 15 1 112-119

(M) INOC IMM 1995 32 NO 8 541-


24/8/95

each 46 bp repeat
containing 30-
recognition site



-889 C-T


IL1α



IL1β (BENSI)


+3953 C-T

-511 C-T




IL1β (CLARK)

Centromere




telomere




IL1β

+8912-9137
5 x 46 bp vnr



IL1α



-889 (C-T)

Primer 1: IL1B - Annealing temperature 55.9°

Primer 2: IL1A - Annealing temperature 53.4°

IL1B

Phalloidin stimulation - -2982 - -2795

LPS stimulation - -3757 - -2729

NFkB sites (within promoter) - -296/-286
- -2761/-2753

29/8/95

Macrophage migration Inhibitory Factor (MIF)

- Exists preformed in macrophages & monocytes
- MIF inhibits TNF secretion
- Inhibits stimulated NO secretion by macrophages.
- 98% identity between mouse and human.
- MIF secretion is induced by TNF α , IFN γ .
- < 1 Kb - 3 exons separated by 2 introns of 189 & 95 bp
- Multiple genes in mouse but not human.

Primers - IL-1 β 3'

Primers were designed using macvector. These primers span 3' AU rich region of IL1 β .

The AU rich region confers instability upon an mRNA sequence - thus cytokines do ~~not~~ have them in order to return cytokine levels to norm as soon as an attack is over (TATTTAT.)

F₁ 5' CAAG CAG AAAA CATG CCC GTC 3' (T_m 57.0)

F₂ 5' GCTCCC ACA TTCT GATG AGCAAC 3' (T_m 57.4)

F₃ 5' CATTC T GATG AGCAA CGCTTC 3' (T_m 56.2)

B₁ 5' TGCAGCACTCAGCAATGAGGAG 3' (T_m 57.8)

30/8/95

Primers spanning the A-T region of IL1 α were also designed

F4 - 5' ATA GCA TAA GTT TCT TGG ACC TCA G 3'

F5/ccc/IL1 β - 5' CAG ATA CTG GAAAA CCA GGC GTAGG 3'

B2/ccc/IL1 α - 5' GCT TGT AGG ACT TGA TTG CAGGTG C 3'

F4/IL1 α /ccc + B2/IL1 α /ccc give a 522 bp fragment

F5/IL1 α /ccc + B2/IL1 α /ccc give a 984 bp fragment

Primers were designed using 'macvector' and checked using BLAST searching to ensure that they matched no other sequence in the human genome

GIBCOBRL

Taq DNA Polymerase

Cat. No. 18038-026

Lot No. **FET401** 500 units; 5 U/ μ l

Exp. Date: **05/97** Store at -20°C (not frost-free).

LICENSED FOR PCR

Description:

Taq DNA Polymerase is isolated from *Thermus aquaticus* YT1. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kD. Taq DNA polymerase is heat-stable and will synthesize DNA at elevated temperatures from single-stranded templates in the presence of a primer. Effective January 1, 1994, Life Technologies™ has modified the unit assay for Taq DNA polymerase, effectively doubling the enzyme concentration for consistent PCR performance.

Components:

18038-026	Taq DNA Polymerase
Y02028	10X PCR Buffer
Y02016	50 mM Magnesium Chloride
90238	1% W-1

Lot No. FET401
Lot No. FD1102
Lot No. FCC102
Lot No. FET403

Unit Definition:

One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

Storage Buffer:

20 mM Tris-HCl (pH 8.0)
0.1 mM EDTA
1 mM DTT
50% (v/v) glycerol
Stabilizers

Unit Assay Conditions:

25 mM TAPS (pH 9.3)
50 mM KCl/
2 mM MgCl₂
1 mM DTT
0.5 mg/ml activated salmon sperm DNA
0.2 mM dATP, dCTP, dGTP, dTTP

10X PCR Buffer:

200 mM Tris-HCl (pH 8.4)
500 mM KCl

The PCR Buffer supplied as a 10X concentrate should be diluted for use.

The 1% solution of the detergent W-1 can be added at a final concentration of 0.05% (v/v) and may improve the thermostability of the enzyme. Store solution at -20°C and thaw at 37°C before use.

Doc. Rev.: 022895

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE IN HUMANS OR ANIMALS.

Printed on Recycled Paper

LIFE TECHNOLOGIES

P.C.R.

dNTP mix .. 10 μ l dATP
10 μ l dCTP
10 μ l dGTP
10 μ l dTTP
360 μ l H₂O \Rightarrow 10mm working stock
400 μ l.

P.C.R. buffer - supplied as 10x \rightarrow working concⁿ = 1x.

1 μ l detergent per 1 μ l Tag.

INITIAL DENATUREⁿ - 96 $^{\circ}$ - 2'

DENATUREⁿ 94 $^{\circ}$ - 1'

56 $^{\circ}$ C Anneal 1'

ELONGⁿ 72 $^{\circ}$ 1'

} X 35

FINAL ELONGⁿ 72 $^{\circ}$ 5'

4 $^{\circ}$

D/S DNA	-	1 O.D.	=	50 μ g/ml
RNA		1 O.D.	=	40 μ g/ml
S/S DNA		1 O.D.	=	33 μ g/ml

P.C.R. 1L1/3 3' END

<u>Reagent</u>	<u>Stock</u>	<u>Use</u>	<u>final</u>
10x P.C.R. buffer	-	5 μ l	1x
MgCl ₂	50mM	2.5 μ l	2.5mM
dNTP's mix	10mM (each)	4 μ l	0.2mM
Taq Pol (GIBCO)	5u/ μ l	0.2 μ l	1u
W-1 (detergent)	-	0.2 μ l	-
Template	50 μ g/ml	2 μ l	100ng/reac ⁿ
Primer mix	2 μ M each	5 μ l	1 μ M
H ₂ O	-	31.1	-
		<u>50μl</u>	

2 drops mineral oil.

F₁ + B₁ annealing Temperature = 56.0°

F₂ + B₁ annealing Temperature = 53.3°

P.C.R. CYCLES

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h53m59s

If a P.C.R. Reaction does not work first time - alter MgCl concⁿ or annealing Temperature or

1/9/95

Primers F1/IL-1~~B~~/CCC
F2/IL1.B/CCC
B1/IL-1~~B~~/CCC

were made and supplied
in dry form. They were
extracted and purified
by ethanol pptⁿ.

MM 737 - B1/IL1.B/CCC
MM 738 - F2/IL1.B/CCC
MM 739 - F1/IL1.B/CCC

	B1	F1	F2
1) Resuspend oligo in 200µl H ₂ O	✓	✓	✓
2) Place in ependorf	✓	✓	✓
3) Add 10µl 3M NaAc.	✓	✓	✓
4) Add 300µl 69% etOH	✓	✓	✓
5) -70°C 1hr. - 11:30am	✓	✓	✓
6) Spin 15' 12,000 rpm	✓	✓	✓
7) Remove Supernatant	✓	✓	✓
8) Wash in 70% ETOH	✓	✓	✓
9) Remove Supernatant	✓	✓	✓
10) Dry in Vacuum ~10'	✓	✓	✓

* When you release vacuum - turn black
mark away from hole - ie: off

Remove hose & turn off H₂O first

* Otherwise you will flood your samples

* AFTER DRYING THE OLIGO SHOULD BE RESUSPENDED IN 500µl *

O.D Programme 10 - Walburg.

step - 10

press RS

Recall press C ⇒

insert blank

H₂O
flush

7

Put tube in touch full press RS

O.D. CALCULATION-

add 170 \rightarrow 200 μ l. - take 5 μ l + 995 H₂O

A₃₂₀ low. .

BI A₂₆₀ 0.3188 = 3

$$= [] = 10.52 \text{ mg/ml}$$

MW = 6787.4

6787.4 g / litre = 1M

6784 mg / ml = 1M

10.52 mg / ml = 1.55 mM

= 1550 μ M \Rightarrow dilute 1:77.5

1 μ l + 77 μ l H₂O

400 μ M solⁿ. \Rightarrow 1 in 4 solⁿ = 195 μ l H₂O + 65 μ l primer

Primer mix - 4 μ l B₁ + 4 μ l F₁
+ 76 152 μ l H₂O

= 260 μ l stock

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
B ₁ 1.0000	0.0191	0.2098	0.3188	0.6362	1.5718	68.957	11.988
F ₁ 2.0000	-0.004	0.1123	0.1435	0.7883	1.2685	68.813	5.0951
P ₂ 3.0000	-0.025	0.1237	0.1154	1.0593	0.9440	124.46	3.4762
H ₂ O 4.0000	-0.020	0.0000	0.0934	0.1770	5.6494	-54.74	6.4125

Column 1

14:13:08 ,30/ 8/95

Run ID : MM 737
 Cycle : 002 UMOL
 End Proc: End CESS (DMT = Off)
 Sequence: Seq01

Average
 Step-wise
 Yield : 98.4
 Total bases = 22
 A= 7, G= 7, C= 5, T= 3, 5= 0, 6= 0, 7= 0, 8= 0
 (mixed bases= 0)

MW: 6787.4

5' > TGC AGC ACT CAG CAA TGA GGA G <3' B1 / 1418 / CCC

Column 2

14:13:09 ,30/ 8/95

Run ID : MM 738
 Cycle : 002 UMOL
 End Proc: End CESS (DMT = Off)
 Sequence: Seq02

Average
 Step-wise
 Yield : 96.8
 Total bases = 23
 A= 6, G= 4, C= 8, T= 5, 5= 0, 6= 0, 7= 0, 8= 0
 (mixed bases= 0)

MW: 6965.6

5' > GCT CCC ACA TTC TGA TGA GCA AC <3' F2 / 1418 / CCC

Column 3

14:13:09 ,30/ 8/95

Run ID : MM 739
 Cycle : 002 UMOL
 End Proc: End CESS (DMT = Off)
 Sequence: Seq03

Average
 Step-wise
 Yield : 96.8
 Total bases = 21
 A= 8, G= 4, C= 7, T= 2, 5= 0, 6= 0, 7= 0, 8= 0
 (mixed bases= 0)

MW: 6390.2

5' > CAA GCA GAA AAC ATG CCC GTC <3' F1 / 1418 / CCC

40 μM

⇒ Stock Solⁿ. (400 μM)

3808.2

1.

546 μM

1.4 solⁿ.91 μL H₂O = 156 μL solⁿ.1:20 dilⁿ.H₂O.

1446

3550.4

2336

3438.4

10770.8

F

$$A_{260} = 0.1435 \times \text{dil factor (1000)}$$

$$= 143.5 \times 33$$

$$= \text{OD of } \underline{\underline{4.735}}$$

$$M_w = 6390$$

$$6390 \text{ g/L} = 1\text{M}$$

$$6390 \text{ mg/mL} = 1\text{M}$$

$$\cancel{6390 \text{ mg/mL} = 1\text{M}}$$

$$4.735 \text{ mg/mL} = 740 \mu\text{M}$$

$$= \underline{\underline{65 \mu\text{L primer} + 130 \mu\text{L H}_2\text{O}}}$$

$$\rightarrow 1:2 \Rightarrow \text{stock sol}^n (400 \mu\text{M})$$

F2

$$A_{260} = 0.1154 \times 1000 =$$

$$115.4 \times 33 = 3808.2$$

$$= 3.808 \text{ mg/mL}$$

$$M_w = 6965$$

$$6965 \text{ g/L} = 1\text{M}$$

$$6965 \text{ mg/mL} = 1\text{M}$$

$$3.808 \text{ mg/mL} = \underline{\underline{546 \mu\text{M}}}$$

$$\Rightarrow \text{to get } 400 \mu\text{M sol}^n \quad 1:1.4 \text{ sol}^n$$

$$\cancel{65 \mu\text{L}} \quad \underline{\underline{65 \mu\text{L primer} + 91 \mu\text{L H}_2\text{O} = 156 \mu\text{L sol}^n}}$$

$$\text{To get from } 400 \rightarrow 20 \mu\text{M} \quad 1:20 \text{ dil}^n$$

$$\text{Primer mix} = 2 \mu\text{L} + 16 \mu\text{L of H}_2\text{O}$$

$$P_1 = B_1 + F_1$$

$$P_2 = B_1 + F_2$$

IL-1B TEST PCR

Two P.C.R reactions were set up 1 with F1/IL1B/ccc + B1/IL1B/ccc (803 bp fragment) and 1 with F2/IL1B/ccc + B1/IL-1B/ccc (403 bp fragment)

Template used = 3 x normal genomic samples

⇒ Set up 8 reactions

Master mix:- Buffer - 40 μ l ✓

MgCl₂ - 20 μ l ✓

dNTP's - 32 μ l ✓

Taq Pol - 1.6 μ l ✓

H₂O - 1.6 μ l ✓

Template

H₂O - 24.88 μ l ✓

344.0 μ l

Pipette 43 μ l of Master mix into each tube ✓

into tubes 1, 2, 3 and 4 add 5 μ l Primer mix 1 ✓

into tubes 5, 6, 7 and 8 add 5 μ l Primer mix 2

tubes 1 and 8 are H₂O controls

lib 0.5:1 Max 0.5 F 0.04 1000 200 of
0.5:1 20 1000 1000 = 1000 1000

17 + 12 = 29

17 + 12 = 29

TUBE No.

*DNA USED

German controls

1/	-	H ₂ O control	+ Primer mix 1	H ₂ O
2/			+ Primer mix 1	Sample No 1
3/			" "	" " 3
4/			" "	" " 7
5/			+ Primer mix 2	(oops Primer " added by mistake)
6/			" "	1
7/			" "	3
8/	-	H ₂ O control	+ Primer mix 2	H ₂ O

Tubes are 'loaded' as above - before adding DNA -

Wax pellets are added to each tube they are then melted @ 75°C for 5' and cooled to 4°C for 3' DNA samples (2 µl) are then added on top of the wax layer before putting on the thermocycler.

No 61 - 75° 5' } wax melting
62 4° 3' }

*The volume of reagent which should be put on top of the wax layer should be 10-30 µl - To ensure equal mixing this can be achieved by adding Master mix to a DNA / H₂O solution which is wax sealed *

IL-12 Primers

DATA

- 394 Synthesis Setup Listing -

(Version 2.00)

Column 1

13:50:56 ,31/ 8/95

Run ID :

MM 740

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq01

Average

Step-wise

Yield : 98.5

Total bases = 25

A= 4, G= 10, C= 3, T= 8, 5= 0, 6= 0, 7= 0, 8= 0
(mixed bases= 0)

MW: 7775.0

11451

5' > GCT TGT AGG ACT TGA TTG CAG GTG G <3'

11422

B2/CCC/IL12

GCT TGT AGG ACTTGA TGG

AGG TGG 3'

ATA GCA TAA GTT TCT GGG ACC TGT

ATA CTG GAA AAC CAG GCG TAG

Column 2

13:50:57 ,31/ 8/95

Run ID :

MM 741

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq02

Average

Step-wise

Yield : 98.8

Total bases = 25

A= 7, G= 6, C= 5, T= 7, 5= 0, 6= 0, 7= 0, 8= 0
(mixed bases= 0)

MW: 7676.0

5' > ATA GCA TAA GTT TCT GGG ACC TCA G <3'

F4/CCC/IL-12

Column 3

13:50:57 ,31/ 8/95

Run ID :

MM 742

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq03

Average

Step-wise

Yield : 98.2

Total bases = 25

A= 9, G= 8, C= 5, T= 3, 5= 0, 6= 0, 7= 0, 8= 0
(mixed bases= 0)

MW: 7742.0

5' > CAG ATA CTG GAA AAC CAG GCG TAG G <3'

F5/CCC/IL1B

52.7° ⇒ 522 bp fragment

54.6° ⇒ 984 bp fragment

11451

529
123

13

IL-12 Primers

DATA

- 394 Auto Analysis Listing -

Time: 16:18:39 , 1/ 9/95

Column #1		Column #2		Column #3	
Seq: Seq01 B2		Seq: Seq02 F4		Seq: Seq03 F5	
Overall: 69.9		Overall: 75.1		Overall: 66.0	
ASWY: 98.5		ASWY: 98.8		ASWY: 98.3	
Num Base	ASWY	Num Base	ASWY	Num Base	ASWY
2 G	100.0	2 A	100.0	2 G	100.0
3 T	97.1	3 C	92.5	3 A	97.3
4 G	94.9	4 T	94.9	4 T	97.7
5 G	96.2	5 C	96.2	5 G	97.6
6 A	96.9	6 C	96.9	6 C	96.2
7 C	95.9	7 A	97.4	7 G	96.8
8 G	96.4	8 G	97.8	8 G	96.7
9 T	96.9	9 G	98.1	9 A	96.6
10 T	97.2	10 G	98.3	10 C	96.9
11 A	96.8	11 T	98.5	11 C	97.2
12 G	97.1	12 C	98.6	12 A	97.5
13 T	97.4	13 T	98.7	13 A	97.7
14 T	97.6	14 T	98.8	14 A	96.8
15 C	97.6	15 T	98.9	15 A	97.1
16 A	97.8	16 G	99.0	16 G	97.3
17 G	97.9	17 A	99.0	17 G	97.4
18 G	98.1	18 A	99.1	18 T	97.6
19 A	98.0	19 T	98.8	19 C	97.7
20 T	98.1	20 A	98.9	20 A	97.8
21 G	98.2	21 C	98.6	21 T	97.9
22 T	98.3	22 G	98.6	22 A	98.0
23 T	98.4	23 A	98.7	23 G	98.1
24 C	98.5	24 T	98.8	24 A	98.2
25 G	98.5	25 A	98.8	25 C	98.3

B2|ccc|142 = 5' GCT TGT AGG ACTTGATGG
AGG TGG 3'

F4|ccc|142 = 5' ATA GCA TAA GTT TCT GGG ACC TGT

F5|ccc|142 = 5' CAG ATA CTG GAA AAC CAGGCG TAG

F4 + B2 - Annealing Temperature = 52.7° ⇒ 522bp fragment

F5 + B2 - Annealing Temperature = 54.6° ⇒ 984bp fragment

11451

529
143

4/9/95

Purify, by ethanol precipitation newly made primers
B2/1412/CCC, F4/1412/CCC, F5/1412/CCC

- 1) Resuspend oligas in 200 μ l Pure H₂O.
- 2) Put 100 μ l into ~~Pete~~ fresh eppendorf
- 3) Add 10 μ l 3M NaAc
- 4) Add 300 μ l EtOH (100%)
- 5) Freeze @ -70°C ~1 hr.
- 6) Spin, 12,000 rpm - 15'
- 7) Remove SN
- 8) Wash pellet in 70% EtOH
- 9) Spin 5' 12,000 rpm
- 10) Dry Under vacuum ~10'
- 11) Resuspend 500 μ l H₂O.
- 12) Measure O.D. (Prog 10. Warburg)

* To put vacuum on - turn black mark away from 'hole'
put on tube and turn on top - to turn off, pull tube
off - then switch off H₂O then lift nozzle *

$$\begin{aligned} B2 \quad A_{260} &= 0.5413 \times \text{dil} \text{ factor}^{(200)} \\ &= 108.26 \times 33 \\ &= 3572.58 \end{aligned}$$

$$\begin{aligned} MW &= 7775 = 7775 \text{ g/l} = 1 \text{ m} \\ &= 7775 \mu\text{g}/\mu\text{l} = 1 \text{ m} \\ 3572.58 &= 460 \mu\text{m} \end{aligned}$$

Dilution factor = 22

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
1.0000	-0.010	0.0104	0.0161	0.7798	1.2824	11.850	0.9111
32 2.0000	0.0630	0.3096	0.5413	0.5156	1.9397	20.491	21.209
F4 3.0000	-0.023	0.1182	0.2347	0.5475	1.8264	23.812	11.122
F5 4.0000	0.0264	0.1240	0.2105	0.5304	1.8853	12.133	8.0636

Run P.C.R. Samples ^(P10111) on 1% Agarose Gel. Add 1 μ l EtBr to Gel.
Make small 30ml gel. 0.5g Agarose
50ml TBE (1x)

- Take 10 μ l PCR product
- Add 4 μ l loading dye

also run 0.5-1 μ l Φ X Hae III marker

Gel 1 loading order.

- 1) Φ X Hae III ~~to~~ Molecular weight marker
- 2) ^{H₂O control} Sample No 1 (Primers B₁/F₁)
- 3) " 3
- 4) " 7
- 5) " 8
- 6)

Gel 2

- 1) Φ X Hae III Molecular wt marker
- 2) ~~No 1~~ ^{H₂O} Sample No 1 (Primers B₁/F₂)
- 3) No 3
- 4) H₂O control

Gel was run @ 80V for ~ 30'.

RESULT

No bands were seen except those of primers - so re-plan expt. Optimise conditions such as magnesium concentration and annealing Temperature. Increasing magnesium concentration lowers the specificity of the reaction as does lowering the annealing temp.

IL1β PCR - MARK II

OPTIMISATION

alter - $MgCl_2$

- Annealing Temp.

$MgCl_2$ - was 2.5 use: 2.0 3.0 3.5

Annealing Temperature was 53' - try 50 47
56' - 53 50

Block 1 Annealing temp: 56°C

Samples 1, 3

Magnesium 2.0 3.0 3.5 μ

Master mix - 1 x 9

Buffer 5 μ (2.5 mM)

$MgCl_2$ 2.5 μ ✓

dNTP's 4 μ ✓

Taq 0.2 μ ✓

W1 0.2 μ ✓

[Template] 2 μ -

Primer 5 μ ✓

H₂O 31.5

50

Mix 2

5 μ (3.5 mM)

3.5 μ

4 μ

0.2 μ

0.2 μ

2 μ

5 μ

30.5

50

Mix 3

5 μ ✓ (4.5 mM)

4.5 μ

4 μ ✓

0.2 μ

0.2 μ

2 μ

5 μ

29

50

$$\begin{array}{r} 31.5 \\ \times 9 \\ \hline 283.5 \end{array}$$

$$\begin{array}{r} 30.5 \\ \times 9 \\ \hline 274.5 \end{array}$$

$$\begin{array}{r} 29.5 \\ \times 9 \\ \hline 265.5 \end{array}$$

	(1) <u>mm1</u>	(2) <u>mm2</u>	(3) <u>mm3</u>
Ⓐ Block 1: -56°	⁽¹⁰⁰⁾⁽¹⁰⁰⁾ 1 2 H ₂ O F₁B₁	1, 2 H ₂ O F₁B₁	1, 2, H ₂ O
Ⓑ Block 2: -53°	1, 2, H ₂ O F₁B₁	1, 2, H ₂ O	1, 2, H ₂ O
Ⓒ Block 3: -50°	1, 2 H ₂ O	1, 2, H ₂ O	1, 2, H ₂ O

O.D. Primers

$$F_4, A_{260} = 0.2347$$

$$0.2347 \times 200 \times 33 = \text{OD}$$

Dilution factor = 10

$$\text{O.D.} = 1549.02$$

$$M_w = 7676.0 =$$

$$1M = 7676 \mu\text{g}/\mu\text{l}$$

$$\Rightarrow 1549 = 202 \mu\text{M}$$

$$F_5, A_{260} = 0.2105 \times 200 \times 33 = \text{OD}$$

$$\text{O.D.} = 1389.3$$

$$M_w = 6390$$

$$= 220 \text{ mM}$$

F₅ - dilution factor = 11
to give 20 mM soln.

$$\boxed{\frac{(A_{260} \times \text{dil}^n \times 0.033 / M_w) \times 10^6}{20} = \text{dil}^n \text{ factor}}$$

$$F_5 \quad \frac{(1.3893) \times 10^6}{20} = 10$$

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~56.0 deg C~~ time: ~~0h 1m 0s~~ slope: ~~2.00 deg/sec~~ time-inc: ~~0 sec~~
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h48m11s

Samples A

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~53.0 deg C~~ time: ~~0h 1m 0s~~ slope: ~~2.00 deg/sec~~ time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h52m59s

Samples B

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~56.0 deg C~~ time: ~~0h 1m 0s~~ slope: ~~2.00 deg/sec~~ time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h59m23s

Samples C

IL-12 PCR**PROVISIONAL.
IL 1a PCR PROTOCOL**

REAGENT	STOCK	USE	FINAL
10 X BUFFER		5 μ l	1X
MgCl ₂	50mM	2.5 μ l	2.5mM
dNTP's	10mM	4 μ l	0.2 mM
Taq Pol	5U/ μ l	0.2 μ l	1U
W-1		0.2 μ l	
TEMPLATE	50 μ g/ml	2 μ l	100ng/react
primer mix	20 μ M each	5 μ l	1 μ M
H ₂ O		31.1 μ l	
		50 μ l	

CYCLES:

INITIAL DENATURATION	96°C	2MINS	
DENATURATION	94 °C	1 MIN	
ANNEALING	F4 (52 °) F5(54 °)	1 MIN	X35
ELONGATION	72°C	1 MIN	
ELONGATION	72°C	5 MINS	
	4°C	INFINITY	

Primer mix 3 = F4/1412/ccc + B2/1412/ccc @ 20mM

Primer mix 4 = F5/1412/ccc + B2/1412/ccc @ 20mM

~~P3~~
~~F4~~
~~B2~~
~~H₂O~~

~~P4~~
~~F5~~
~~B2~~
~~H₂O~~

5/9/95

Run P.C.R. samples on a 1% ^{Agarose} ~~Agarose~~ Gel

10 μ l P.C.R. products

4 μ l bromophenol blue

load 14 μ l

Loading order

- | | | | |
|-----|--------------------------------------|-----|-------------------------------------|
| 1) | 1 μ l Φ X Hae III Mw marker | 18) | 1 μ l Φ X Hae III |
| 2) | 14 μ l A ₁ (1) | 19) | 14 μ l B ₃ (1) |
| 3) | " A ₁ (2) | 20) | " B ₃ (2) |
| 4) | A ₁ H ₂ O | 21) | B ₃ (H ₂ O) |
| 5) | A ₂ (1) | 22) | " C ₁ (1) |
| 6) | A ₂ (2) | 23) | " C ₁ (2) |
| 7) | A ₂ (H ₂ O) | 24) | " C ₁ (3) |
| 8) | A ₃ (1) | 25) | " C ₂ (1) |
| 9) | A ₃ (2) | 26) | " C ₂ (2) |
| 10) | A ₃ (H ₂ O) | 27) | " C ₂ (H ₂ O) |
| 11) | B ₁ (1) | 28) | " C ₃ (1) |
| 12) | B ₁ (2) | 29) | " C ₃ (2) |
| 13) | B ₁ (H ₂ O) | 30) | " C ₃ (H ₂ O) |
| 14) | B ₂ (1) | 31) | " Φ X Hae III |
| 15) | B ₂ (2) | | |
| 16) | B ₂ (H ₂ O) | | |
| 17) | 1 μ l Φ X Hae III Mw marker | | |

Run gel at 120 V. for ~ 1/2 hour

Set up 16 12 PCR as per provisional protocol P19.

<u>mm(1)</u>			<u>mm(2)</u>		
buffer	-	15 μ l	buffer	-	15 μ l
MgCl ₂	-	7.5 μ l	MgCl ₂	-	7.5 μ l
dNTP's	-	12 μ l	dNTP's	-	12 μ l
Taq Pol	-	0.6 μ l	Taq	-	0.6 μ l
W-1	-	0.6 μ l	W1	-	0.6 μ l
Primer mix	(F4/B2)	6 μ l	Primer mix	(F5/B2)	6 μ l
H ₂ O		93.3 μ l	H ₂ O		93.3 μ l
		150 μ l			150 μ l

For each primer pair use two different DNA samples
and 1 water control

FURTHER OPTIMISATION OF IL1B PCR

MgCl₂ - 1.0 μl 1.5 μl 5.0 μl
Temp. 57 °C 60 °C 63 °C

mm ₁	mm ₂	mm ₃
Extr - 5 ✓ (45)	5 μl (45)	5 μl
MgCl ₂ - 2.5 1.0 (9)	1.5 μl (13.5)	2.0 μl
dNTPs - 4 ✓ (36)	4 μl (36)	4 μl
Taq - 0.2 (1.8)	0.2 μl (1.8)	0.2 μl
W ₁ - 0.2 (1.8)	0.2 μl (1.8)	0.2 μl
Primer - 5 (45)	5 μl (45)	5 μl
H ₂ O - 32.6 (23.4)	32.1 μl (288.9)	288.6 μl (284.4)
<u>48 μl</u>	<u>48</u>	<u>48 μl</u>

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~65.0 deg C~~ 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h47m59s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~65.0 deg C~~ 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: ~~65.0 deg C~~ 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h38m59s



Loading order.

2)	A ₁ (1)	Temp 57°C	mm 1	DNA sample 96
3)	A ₁ (2)	"	mm 1	" 97
4)	A ₁ (H ₂ O)	"	mm 1	H ₂ O control
5)	A ₂ (1)	"	mm 2	" 96
6)	A ₂ (2)	"	↓	97
7)	A ₂ (H ₂ O)	"	↓	H ₂ O control
8)	A ₃ (1)	"	mm 3	96
9)	A ₃ (2)	"	↓	97
10)	A ₃ (H ₂ O)	"	↓	H ₂ O control

11)	B ₁ (1)	Temp 60°	mm 1	# 96
12)	B ₁ (2)	"	↓	# 97
13)	B ₁ (H ₂ O)	"	↓	H ₂ O control
14)	B ₂ (1)	"	mm 2	# 96
15)	B ₂ (2)	" X	↓	# 97
16)	B ₂ (H ₂ O)	"	↓	H ₂ O control
17)	B ₃ (1)	"	mm 3	# 96
18)	B ₃ (2)	"	↓	# 97
19)	B ₃ (H ₂ O)	"	↓	H ₂ O control



GA)

20)	C ₁ (1)	Temp 63°	mm 1	# 96
21)	C ₁ (2)	"	↓	# 97
22)	C ₁ (H ₂ O)	"	↓	H ₂ O control
23)	C ₂ (1)	"	mm 2	# 96
24)	C ₂ (2)	" X	↓	# 97
25)	C ₂ (H ₂ O)	"	↓	H ₂ O control
26)	C ₃ (1)	"	mm 3	# 96
27)	C ₃ (2)	"	↓	# 97
	C ₃ (H ₂ O)	"	↓	H ₂ O control

* ALL PCR REACTIONS ARE BEING CARRIED OUT IN 25µL TO MINIMISE WASTAGE OF REAGENTS *

To ensure that all my reagents are working and I am not doing anything silly I will try out an already proven P.C.R. - (Alisons)

Date	5/9/95.
Number of Samples	3
Disease	Any!
P'morph	Taq 1
Water	48.8 μ l
Buffer	8 μ l
Magnesium	8 μ l
dNTPs	6.4 μ l
Primers	8 μ l
Taq	0.4 μ l
W-1	0.4 μ l
Template	2 μ l



1b1 B RESULTS Gel loading order P23 -

A number of bands, of the correct size were seen. Those most prominent were:

- 60°C Annealing - 1.5 μ l $MgCl_2^{2+}$
 - * - 63°C Annealing - 1.5 μ l $MgCl_2^{2+}$
- For Gel pictures see Above

To ensure that
I am not doing
already proven

ing are
ut an

Date 5./9./95.

Number of
Samples 3

Disease Any!

P'morph Taq 1

Water 48.8 μ l

Buffer 8 μ l

Magnesium 8 μ l



dNTPs 6.4 μ l

Primers 8 μ l

Taq 0.4 μ l

W-1 0.4 μ l

Template 2 μ l



1.1 B RESULTS Gel loading order P23

A number of bands, of the correct size were seen.
Those most prominent were:

- 60°C Annealing - 1.5 μ l $MgCl_2^{2+}$
- * - 63°C Annealing - 1.5 μ l $MgCl_2^{2+}$

For Gel pictures see Above

OPTIMISATION OF IL1 α PCR

Primers F4/B2

	mm1	mm2	mm3	mm4	mm5	mm6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Mg	0.5 ^{4.5}	1 ⁹	1.5 ^{13.5}	2 ¹⁸	2.5 ^{22.5}	3 ²⁷
dNTPs	2	2	2	2	2	2 (15)
Primers	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Taq	0.1 μ l	0.1	0.1	0.1	0.1	0.1 (0.9)
W-1	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
H ₂ O	17.3 ^{155.7}	16.8 ^{151.2}	16.3 ^{146.7}	15.8 ^{142.2}	15.3 ^{137.7}	14.8 ^{133.2}
	25	25	25	25	25	25

Predicted annealing Temp: - 53°C

DNA (Psoraleis) sample: X6 + H₂O

X6 + H₂O

H₂O

Make 3x amount - for studies at diff Temperatures
ie a total of 9x above

buffer.

Carry out above PCR at same MgCl₂ concn
but at 50 and 56°C

53°C set - A₁-6 + A H₂O } labelled 50°C on side
B₁-6 + B H₂O }

56° set: A₁-6 + A H₂O } labelled 56°C on side
B₁-6 + B H₂O }

Summary of ILI-2 OPTIMISATION

6 different master-mixes @ 6 diff Mg concn.

1) 0.5 μ l

4) 2.0 μ l

2) 1 μ l

5) 2.5 μ l

3) 1.5 μ l

6) 3.0 μ l

⇒ Each master mix was tested out with 2 DNA samples and water controls

⇒ Reactions were carried out at 50, 53 and 56°C

Loading order IL-12 OPTIMISATION

1) M.w Marker Q

2)	DNA sample	125	mm 1	50°C
3)	"	135	2 1	"
4)		125	3 2	"
5)		135	4 2	"
6)		125	5 3	"
7)		135	3	"
8)		125	4	"
9)		135	4	"
10)		125	5	"
11)		135	5	"
12)		125	6	"
13)		135	6	"
14)	H ₂ O control		3	"
15)	H ₂ O control		3	"



1) M.w Marker.

2)	Sample	125	mm 1	53°C
3)		135	1	"
4)		125	2	"
5)		135	2	"
6)		125	3	"
7)		135	3	"
8)		125	4 - 2005	"
9)		135	4	"
10)		125	5	"
11)		135	5	"
12)		125	6	"
13)		135	6	"
14)		125	3	"
15)		135	6	"

1) Molecular wt marker

2)	ONA sample	125	mm 1	56°C
3)	"	135	" 1	"
4)	"	125	" 2	"
5)	"	135	" 2	"
6)	"	125	" 3	"
7)	"	135	" 3	"
8)	"	125	" 4	"
9)	"	135	" 4	"
10)	"	125	" 5	"
11)	"	135	" 5	"
12)	"	125	" 6	"
13)	"	135	" 6	"
14)	"	125	"	"
15)	"	135	"	"
14)	H ₂ O control		" 6	"
15)	H ₂ O control		" 3	"

RESULTS




at 50°C only smears were
observed

Conditions for ILI & PCR are high Temp
Low magnesium
53°C high magnesium

1) Molecular wt marker

2)	ONA sample	125	mm 1	56°C
3)	"	135	" 1	"
4)	"	125	" 2	"
5)	"	135	" 2	"
6)	"	125	" 3	"
7)	"	135	" 3	"
8)	"	125	" 4	"
9)	"	135	" 4	"
10)	"	125		
11)	"	135		
12)	"	125		
13)	"	135		
14)	/	/	"	12
15)	/	/	"	13
16)	H ₂ O control			
17)	H ₂ O control			

RESULTS



at 50°C only smears were observed

From Above Gel it seems that best conditions for PCR are high Temp low magnesium
53°C 1mM magnesium

6/9/95

OPTIMISATION OF IL1 β F2/B1 PRIMERS

Tried already:- 53°C annealing Temp.
2.5mM Magnesium.

Try now Mg 1 2 3 4 5 6 mM
Temp 50° 56° 60°

Reagent	1	2	3	4	5	6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Mg	0.5	1	1.5	2	2.5	3 (27)
dNTPs	2	2	2	2	2	2 (18)
Primers	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Taq	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
W1	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
H ₂ O	16.3 ^(146.7)	15.8 ^(142.2)	15.3 ^(137.7)	14.8 ^(135.2)	14.3 ^(129.7)	13.8 ^(124.4)
	24	24	24	24	24	

label samples:-

DNA sample - 157 = 7
- 167 = 4

IL-1B OPTIMISATION Primers F2/B1

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 50.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 3h 2m23s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h52m11s

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 60.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h46m11s

Loading Order

	<u>DNA sample</u>	<u>Master mix</u>	<u>Temperature</u>
①	Mw marker		
②	157	mm 1	50°
③	167	1	50°
④	157	2	"
⑤	167	2	"
⑥	157	3	"
⑦	167	3	"
⑧	157	4	"
⑨	167	4	"
⑩	157	5	"
⑪	167	5	"
⑫	157	6	"
⑬	167	6	"
⑭	157 H ₂ O		"
⑮	167 H ₂ O		"


①	157		
②	157	mm 1	56°
③	167	1	"
④	157	2	"
⑤	167	2	"
⑥	157	3	"
⑦	167	3	"
⑧	157	4	"
⑨	167	4	"
⑩	157	5	"
⑪	167	5	"
⑫	157	6	"

	<u>DNA sample</u>	<u>Master Mix</u>	<u>Temperature</u>
⑬	167	6	56°
⑭	157 H ₂ O	6	"
⑮	H ₂ O	3	"
⑯			

①			
②	157 (X)	mm ₁	60°
③	167 (X)	1	60°
④	157	2	60°
⑤	167	2	60°
⑥	157	3	60°
⑦	167	3	60°
⑧	157	4	60°
⑨	167	4	60°
⑩	157	5	60°
⑪	157	5	60°
⑫	156	6	60°
⑬	167	6	60°
⑭	H ₂ O	6	60°
⑮	H ₂ O	3	60°
⑯			

- 1st gel didn't run properly - Will try to Re-run with remaining samples. - 56° low Mg looked Good.
 ⇒ Predicted Annealing Temp = 53°.
 Re-run gel with remainder of samples ⇒ See over.

RESULTS



IL1B F2/B1 Primers

OPTIMISED CONDITIONS FOR IL-1B (F2/B1)
 = 2mM magnesium
 56°C annealing Temp.

loading order: -

			Temp.
1)	X174 Hae III	MW marker.	
2)	157 (Francis Borais)	mm1	56°
3)	167	mm1	"
4)	157	" 2	"
5)	167	2	"
6)	157	3	"
7)	167	3	"
8)	157	4	"
9)	167	4	"
10)	157	5	"
11)	167	5	"
12)	157	6	"
13)	167	6	"

IL-1B FURTHER OPTIMISATION. (F./B.)

⇒ Previously shown that 63°C at 1.0-1.5 mM magnesium produced best results. Will try.

Temperatures :- 62°/64°/66°

MgCl₂

0.5/1/1.5

1/1.5/2

	MASTER-MIX 1	MASTER MIX 2	MASTER MIX 3
BUFFER	2.5µl	2.5µl	2.5µl (15)
Mg	0.25 ⁽²⁵⁾ ₍₃₎	0.5 ⁽¹⁵⁾ _{4.5}	0.75 ⁽⁶⁾
dNTP's	24µl ¹⁰	24µl	24µl (12)
primer	2.5µl	2.5µl	2.5µl (15)
Taq	0.1	0.1	0.1 (0.6)
w-1	0.1	0.1	0.1 (0.6)
H ₂ O	18.8 ^(12.9)	16.05 ^(16.3)	15.8 ⁽¹⁹⁾ 14.8
TOTAL	24µl	24µl	24µl

DNA samples used Psoriasis 164 - P

174 - Q

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 66.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h37m47s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 64.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h40m11s

LOADING ORDER

1) QX174 HaeIII

2)	Psoriasis	sample No.	164	mm1	62°
3)	"	"	174	1	62°
4)	"	"	164	2	62°
5)	"	"	174	2	62°
6)	"	"	164	3	62°
7)	"	"	174	3	62°
8)	"	H ₂ O control		1	62°
9)	"	"	164	1	64°
10)	"	"	174	1	64°
11)	"	"	164	2	64°
12)	"	"	174	2	64°
13)	"	"	164	3	64°
14)	"	"	174	3	64°
15)	"	H ₂ O control		2	64°
16)	"	"	164	1	66°
17)	"	"	174	1	66°
18)	"	"	(164	2	66°
19)	"	"	174	2	66°
20)	"	"	164	3	66°
21)	"	"	174	3	66°
22)	"	H ₂ O control		3	66°
23)	QX174	HaeIII			

From Yesterdays disastrous Cal:-


1) Q x 174 Hae III

2)	157	1	56°
3)	167	1	56°
4)	157	2	56°
5)	167	2	56°
6)	157	3	56°
7)	167	3	56°
8)	157	4	56°
9)	167	4	56°
10)	157	5	56°
11)	167	5	56°
12)	157	6	56°
13)	167	6	56°
14)	157	1	60°

RESULTS

603

310




IL1B F2/B1 Primers

Best band = 2mM mg
56° Annealing Temp

03

603



IL1B F1/B1

Best band = 66°C

MgCl²⁺ = 2mM

wrong sized fragment
must have used
wrong Primers

7/8/95

IL-12 - F5/B2 OPTIMISATION

Predicted annealing Temperature 54.6°

	Mix 1	Mix 2	Mix 3	Mix 4	Mix 5	Mix 6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5
MgCl ₂	0.5	1	1.5	2	2.5	3.0
dNTP's	2	2	2	2	2	2
Primers	2.5	2.5	2.5	2.5	2.5	2.5
Taq	0.1	0.1	0.1	0.1	0.1	0.1
W-1	0.1	0.1	0.1	0.1	0.1	0.1
H ₂ O	14.8	16.05	15.3	14.8	14.8	13.8
	24µl	24µl	24µl	24µl	24µl	24µl

Try above mixes at 52°C / 54°C / 56°C

Make up enough for 8 reactions

	Mix 1	Mix 2	Mix 3	Mix 4	Mix 5	Mix 6
Buffer	20	20	20	20	20	20
MgCl ₂	4	8	12	16	20	24
dNTP's	16	16	16	16	16	16
Primers	20	20	20	20	20	20
Taq	0.8	0.8	0.8	0.8	0.8	0.8
W-1	0.8	0.8	0.8	0.8	0.8	0.8
H ₂ O	130.4	128.4	122.4	118.4	114.4	110.4
	192	192	192	192	192	192

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 52.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h58m47s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
~~step 2: temp: 54.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec~~
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h55m47s

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
~~step 2: temp: 54.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec~~
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h52m11s

loading Order

1) 2X1.74 marker

2) 157

1

56°

3) 167

1

56°

4) 157

2

56°

5) 167

2

56°

6) 157

3

56°

7) 167

3

56°

8) 157

4

56°

9) 167

4

56°

10) 157

5

56°

11) 167

5

56°

12) 157

6

56°

13) 167

6

56°

14)

QX174 Hae II Molecular wt marker

1)	157	mm1	52°
2)	167	mm1	52°
3)	157	mm2	52°
4)	167	mm2	52°
5)	157	mm3	52°
6)	167	mm3	52°
7)	167	mm4	52°
8)	167	mm4	52°
9)	157	mm5	52°
10)	167	mm5	52°
11)	157	mm6	52°
12)	167	mm6	52°

111

112


113

114

157	mm1	54°
167	mm1	54°
157	mm2	54°
167	mm2	54°
157	mm3	54°
167	mm3	54°
157	mm4	54°
167	mm4	54°
157	mm5	54°
167	mm5	54°
157	mm6	54°
167	mm6	54°



Results.

54°



From these results
the best bands
are appearing at low
MgCl₂ and 54°C.

⇒ Try 1mM mg
53 & 55°C
+ 2µl DNA sample.



← DISASTER
GEL!!!

IL-1B Fl / B1 OPTIMISATION

Temp 62° 64° 66°

Mg 1 1.5 2.


Reagent	mm1	mm2	mm3
Buffer	2.5	2.5	2.5 (20)
Mg	0.5 (4)	0.75 (6)	1 (8)
dNTPS	2	2	2 (16)
Primers	2.5	2.5	2.5 (20)
WT	0.1	0.1	0.1 (0.8)
Taq	0.1	0.1	0.1 (0.8)
H ₂ O	18.8 (112.8)	16.05 (96.3)	15.8 (94.8)
	24	24	24

Use 2µl DNA Template

Make up enough for 8 reactions (quantities in brackets)

loading Order.

1)	Φ X174 Hae III		62
2)	157	mm 1	62°
3)	167	mm 1	62°
4)	157	mm 2	62°
5)	167	mm 2	62°
6)	157	mm 3	62°
7)	167	mm 3	62°
8)	157	mm 2	64°
9)	167	mm 2	64°
10)	157	mm 3	66°
11)	167	mm 3	66°
12)	^{H1 H2 H3} ΦX174 Mu2 marker		



B₁/F₁

* best sample = 66°C Annealing
2mm Mg

- 803 bp.

It seems that for the larger PCR products -
using 2μl DNA works best.

8/9/95 1L-12 P.C.R. OPTIMISATION.

F5/B2 (P38 For initial expts)

Try 2mm magnesium at 58°, 55° and 60°

Reagents	Per reaction	Master mix (12)	tubes labelled
- Buffer	2.5	30	
- MgCl ₂	1.0	24	Sp) R, H) S,
- dNTPs	2	24	
- Primers	2.5	30	
- w1	0.1	1.2	
- Taq	0.1	1.2	
- H ₂ O	14.8	177.6	
TOTAL	24.0		

9d 200

F4/B2 (INITIAL OPTIMISATION P27)

First results not convincing: - Try

1mM 2mM 3mM MgCl₂
Annealing Temp 58° 60° 62°.

REAGENTS	mm1	mm2	mm3
Buffer	2.5	2.5	2.5 (20)
MgCl ₂	0.5 (4)	1.0 (8)	1.5 (12)
dNTPs	2	2	2 (16)
Primers	2.5	2.5	2.5 (20)
w-i	0.1	0.1	0.1 (0.8)
Taq	0.1	0.1	0.1 (0.8)
H ₂ O	16.3 (120.4)	15.8 (126.6)	15.3 (122.4)
TOTAL	24.0	24.0	24.0

make enough for 8 reactions

L = 157

m = 157

14/11/95


2x F1 PCR's using biotinylated primers (for Mohammed)

Boff	
buffer	210 μ l
MgCl ₂	84 μ l
dNTP	168 μ l
Primer	F1 105 R1 21
W1	4.2
Taq	4.2
H ₂ O	1345 μ l

Do 2 Samples 140/188 with paraffin oil
- 2 Samples " " with Wax

Running order: -

- 1) QX HaeIII
- 2) 140 oil
- 3) 188 oil
- 4) 140 wax
- 5) 188 wax



20/11/95

⇒ PCR up a further (5x10g) reactions for two samples. This time to avoid any ethanol ppt steps the sample will be run straight through a column after band purification - this will miss out any ethanol & also remove excess primer

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 35 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 66.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 35 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 2h37m47s

Prepare 2 samples with paraffin oil and two with wax, incase paraffin oil is interfering with the reaction

28/11/95

SSCP - Try IL1 α -889 PCR cold to see if it works. - if it does redo hot & follow Taras protocol

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 45 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 45 times
step 4: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 6: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 3h36m 8s

Date...../...../..... 58°C

Number of Samples 10

Disease.....

P'morph.....

Water	6.1 μ l
Buffer	1.0 μ l
Magnesium	1.0 μ l
dNTPs	0.8 μ l
Primers M.x	1.0 μ l
Taq	0.05 μ l
W-1	0.05 μ l
Template μ l

IL1 α -889 (30 μ l PCR)

each tube:

3 μ l 10x Buffer

2.4 μ l dNTPs

1.5 μ l 603 primer

1.5 μ l 604 primer

1.2 μ l MgCl₂

0.15 μ l W-1

0.15 μ l ~~W-1~~ Taq

18.9 μ l H₂O

28.8 μ l

add 28.8 μ l to each tube + 1.2 μ l DNA

30/11/95

Try PCR using Angie / Nicolas Protocol. - and
annealing temperature at 50°C

Make up two times the volume - 5 DNA
Samples + 1 H₂O control.

10 µl PCR

100

1 µl buffer.

0.8 µl dNTP

0.5 µl each Primer $\begin{matrix} 603 \\ 604 \end{matrix}$

0.4 µl MgCl₂

0.05 µl w1

0.05 µl Taq

6.3 µl H₂O.

8/1

1996

Biotinylated primers are correctly biotinylated - START AGAIN

Re-do PCR using B primer only

When doing PCR do 45 cycles instead of 35
this will also help insure against excess primers
interfering with the reaction.

Do 4 PCR'S (5 x 100 μ l reactions) 2 normal
and two out out of gel & purifield

REAGENT	VOLUME	Samples used
Buffer	140 200	
dNTP	88 160	133
primer	55 μ l each 100 μ l ea.	139
MgCl ₂	44 80	173
Taq	44 8	179
W-1	44 8	
H ₂ O	695.2 μ l 1264	
	1100 μ l	

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 60.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h28m21s

yes PCR

Primer	50 each	- PCR O.D. NOT
MgCl ₂	40	WORK - TRY
Taq	5	AGAIN WITH
W-1	5	FRESH REAGENTS
H ₂ O	630 μ l	esp. dNTP

2/1/96

Re-dissolve ^{returned} primers in 150 μ l H₂O - measure O.D.

5 μ l \rightarrow 995 μ l H₂O \rightarrow O.D.

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
1.0000	0.0036	0.1146	0.1865	0.6065	1.6488	33.653	7.5100
2.0000	0.0048	0.1625	0.2776	0.5782	1.7296	38.189	11.478

B1 - dilution factor =

$$A_{260} \times \text{DILN} \times 0.033 / \text{WM} \times 10^6$$

$$= \left(\frac{0.1865 \times 200 \times 0.033}{67876} \right) 10^6$$

20

1:9 dilution

B2

$$\left(\frac{0.2776 \times 200 \times 0.033}{7775} \right) 10^6$$

20

= 1:9 dilution

1:11 dilution

Do ILIP/PCR's

	M1	M2
Buffer	200	200
MgCl ₂	80	80
dNTPs	160	160
Primer	100ul $\begin{matrix} F2 \\ B1 \end{matrix}$	100ul $\begin{matrix} F5 \\ B2 \end{matrix}$
Taq	8	8
W-1	8	8
H ₂ O	1264	1264
	2,000	2,000

4 reactions @ 100ul x 5

Use samples $\begin{matrix} 115 \\ 151 \end{matrix} \} 1,1 \text{ for } +3953$

$\begin{matrix} 121 \\ 126 \end{matrix} \} 2,2 \text{ for } +3953$


23/1/96

Running order:-

Run 10ul on

1% Agarose gel

- 1) QX Hae III
- 2) 115 F2
- 3) 151 F2
- 4) 121 F2
- 5) 126 F2
- 6) 115 F5
- 7) 151 F5
- 8) 121 F5
- 9) 126 F5
- 10) ~~QX~~ H₂O control
- 11) QX Hae III



Do 115 and 121 for each Normally
- 151 and 126 extract band from gel to remove primers

RUNNING ORDER - (ON TAE GEL)

F5 151

F2 151


F2 126

F5 126

- Bands were cut out of the gel, and purified through glass wool and ethanol precipitation

Samples were made single-stranded using dynabeads. Samples of supernatants were kept & run on a 1% Agarose gel.

RUNNING ORDER



- 1) QX Hae III
- 2) 1st wash F2 121
- 3) " " F5 151
- 4) 1st wash SS DNA F2 121
- 5) SS DNA 121
- 6) 1st wash F2 121
- 7) " " 115
- 8) SS DNA F5 151

We can see from the previous gel picture that for samples which had excess primer removed by running samples down the gel showed no loss of sample on 1st wash and a clear band appeared on the lane incorporating ss DNA. - The band may not be very bright but it is well known that ss DNA does not incorporate ss DNA so well.


⇒ Try sequencing F5 151 DNA.

23/1/96

Sequencing of PCR product (not biotin)
 is now possible and very successful.
 prepare PCR's for 10 samples
 F₁ - B₁ & F₅ - B₂ 5 Tag 1,1
 and 5, x Tag 2,2.

⇒ enough for 100 reactions each

⇒	Reagent	MM1	MM2
	Buffer	500µl	500µl
	MgCl ₂	200µl ✓	200µl
	dNTP	400µl	400µl
	Primer	250µl each	250µl each
	w-1	20µl	20µl
	Taq	20µl	20µl
	H ₂ O	2960	2960
		4600	4600



Running order

- 1) 115
- 2) 121 ✓
- 3) 126 ✓
- 4) 144 ✓
- 5) 165 ✓
- 6) 113
- 7) 139 ✓
- 8) 157 ✓
- 9) 112
- 10) 157 ✓
- 11) 120

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc:

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc:

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc:


end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0

total runtime (approx.): 3h36m 1s

Take all samples that worked, pool their DNA and ethanol extract - then run on a gel and band purify the sample for sequencing



← Run on a gel to check samples had not been lost by above process


1/2/95

MUTS



0/

- 1) Take 50 μ l DNA (genomic) [LIX 50] 96 $^{\circ}$ C - 5'
- 2) Add 140 μ l buffer ^{TRIS HCL PH 7.5} PCR BUFFER - to make volume 190 μ l 37 $^{\circ}$ C - 5'
- 3) Add 10 μ l beads X3
- 4) Agitate gently by flicking - incubate @ 37 $^{\circ}$ C - 1/2 hr with shaking
- 5) Centrifuge briefly
- 6) Apply magnet
- 7) Wash in 200 μ l Wash buffer X2
- 8) add buffer from step 2 - heat - 75 $^{\circ}$ C - 15' Remove supernatant - & PCR - Using FI/BI primers as before! - Run on gel to see if anything's there!



1/2/96

- Do a 100 μ l PCR reaction with Alisons Tag primers expected fragment size

- also Do 100 μ l PCR with IL1 β F2-B, expected fragment size - 400 bp

(Tag)	Reagent	20 μ l	
		Volume	MM
	Buffer	12.2 2	20
	MgCl ₂	2	20
	dNTP	2 1.6	16
	WT	0.1	1
	Tag	0.1	1
	Primer	2 μ l	20
	H ₂ O	12.2	122

1.6
x 12

3.12
160

198

(IL1 β)	Reagent	Volume
	Buffer	15 μ l 20 μ l
	MgCl ₂	6 μ l 8 μ l
	dNTP	16 μ l
	Primer	10 μ l each
	Tag	0.8 μ l
	WT	0.8 μ l
	H ₂ O	120 μ l

10 μ l template

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
begin cycle 45 times:
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec
end cycle 45 times
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec
total runtime (approx.): 3h 36m 1s

5/2/96

Samples for sequencing:-

1 - 121

2 - 126

3 - 157

4 - 151

5 - 144

6 - 139

7 - 165

- Give Hazel F2 primer to sequence
half of each of the samples:- Tag FS sequencing

Do PCR's as before ←

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h36m 1s

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 64.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h24m 27s

90

Do PCR's on IL1 β for direct sequencing from PCR product.


Do PCR with F5/B2 and sequencing using both F5/F4 primers

F1/B1 PCR - again both F1/F2 primers

	REAGENT	Volume	Volume
10 reactions x 10	buffer	500 μ l	500
= 100	dNTP's	200 400	400
	MgCl ₂	200	200
	primer	250 each	250 each
	Tag	20 μ l	20 μ l
	W-1	20 μ l	20 μ l
	H ₂ O	2960 μ l	2960 μ l
		4600 μ l	4600 μ l

Running order:-

- 1) QX Hae III
- 2) 115 F⁻
- 3) 151
- 4) 157
- 5) 165
- 6) 174
- 7) 121
- 8) 126
- 9) 189
- 10) 144
- 11) 152
- 12)




3/2/96

- Run PCR products on gel - band extract and
ethanol precipitate - Give to Hazel to
sequence

Gel to check products are still there

Running Order

- 1) QX Mac III marker
- 2) F5 174
- 3) F1 126
- 4) F5 151
- 5) F1 151
- 6) F1 174
- 7) F5 152
- 8) F5 157
- 9) F1 139
- 10) F5 121
- 11) F5 139
- 12) F5 165
- 13) F1 + 152
- 14) F1 165
- 15) F5 115
- 16) F5 126
- 17) F1 157
- 18) F1 144
- 19) F1 115




10/2/96

Do Biotinylated PCR's for 5x 1,1
5x 2,2.

REAGENT	VOLUME
buffer	500 μ l
MgCl ₂	200 μ l
dNTP	400 μ l
Primer	250 μ l each (Used 200 F1 100 BB1)
W-1	20 μ l
Taq	20 μ l
H ₂ O	2960 μ l

Run samples on a 1% agarose gel



RUNNING ORDER

- 1) QX Hing
- 2) 115 F1/B61
- 3) 151 F1/B61
- 4) 157 F1/B61
- 5) 165 F1/B61
- 6) 174 F1/B61
- 7) 121 F1/B61
- 8) 126 F1/B61
- 9) 139 F1/B61
- 10) 144 F1/B61
- 11) 152 F1/B61


11/12/06

Pool all 5 x 100 μ l's and ethanol
precipitate resuspend in 20 μ l for
running ~~off~~ on TAE gel & band extract

Running Order

- 1) 139
- 2) 151
- 3) 126
- 4) 144
- 5) 115
- 6) 121
- 7) 157
- 8) 174
- 9) 152
- 10) 165

After band extraction, samples were
ethanol precipitated and resuspended in 40 μ l
H₂O. - 1 μ l of sample was run on a 1%
agarose gel to check that I have not lost
them during ethanol precipitation.



14/2/96

Do IL-1 Tag PCR on 50 patients (4x100ul.)
for control experiments using muts 5 beads
⇒ Worked OK.

Mohammeds Solutions -FILS7
-FI-174
-FI-152
-FI-121
-FI-139

$$1\% = 1\text{ml} / 100$$

$$0.005\% =$$

$$0.005\%$$

$$5\mu\text{l} \rightarrow 100\text{ml}$$

$$0.5 \rightarrow 10\text{ml}$$

$$20\mu\text{l} \rightarrow 10\text{ml}$$

$$1\mu\text{l} \rightarrow 40$$

15/2/96

MUTS - New Protocol

- Carry out reaction in 10mM Tris-HCl (Not IM)
- Put ^{0.005%} tween 20 in buffer rather than BSA.
- ensure buffer is FILTERED through 0.45 μ m
- Binding occurs at room temp. Not 37°C
- Wash beads with 10mM T/HCl (Not tween)

CONTROLS

- (1) Beads - No DNA - boil, remove 5' - PCR as a H₂O control.
 - (2) One @ 37°C and one at Room temp
- beads with no bsa - but tween were OK

REACTION BUFFER

- 1M Tris HCl pH 7.5
- 250mM MgCl₂
- 5mM dTT
- 0.5mM EDTA

50X

← PCR buffer!

+ add MgCl₂

~~1/10~~ Equal Volume

Block beads with 5mM P.V.P.


2g → 10mls

15/2/96

- Again samples were heated to 75°C for fifteen minutes and supernatant removed immediately

21/2/96


Sequencing has revealed a possible G-C Change.
 Design method for screening possible polymorphs.
 - No enzyme (except Tsp RI) cuts at the altered site. Therefore a cut site must be engineered - one possibility is to engineer a cut site for Hinf! in a primer. There are 2 intrinsic Hinf sites + the engineered one.



EXPECTED BANDING PATTERNS

	<u>GG</u>	<u>CC</u>	<u>GC</u>
89	—		— ? faint
76	—	—	—
61	—	—	—
54		—	—
35		—	—) faint

Design oligo which creates the hinf site.
~~Reverse~~



21/ 2/96

Column #3

MM828

Seq: Seq03

Overall: 50.5

ASWY: 98.0

Num Base ASWY

2 C	100.0
3 C	98.5
4 A	96.0
5 T	97.0
6 T	97.6
7 T	98.0
8 A	98.0
9 A	98.0
10 A	97.9
11 T	98.0
12 A	98.1
13 A	97.4
14 T	97.6
15 C	97.8
16 T	97.9
17 G	98.0
18 A	98.1
19 G	98.2
20 C	98.3
21 T	98.2
22 T	97.4
23 A	97.5
24 T	97.7
25 A	97.8
26 T	97.8
27 T	97.9
28 C	98.0
29 A	98.1
30 T	98.1
31 T	97.7
32 G	97.8
33 A	97.9
34 G	98.0
35 T	98.0

SEQUENCE +6911 Primer

5' TGA GTT TTA TAT ATT ATT
CGA GTC TAA TAA ATT TACC 3'

Predicted annealing temperature
57°C

(Using hypercard)

Column 1

14:53:21 ,23/ 2/96

Run ID : NM831
Cycle : 002 UMOL
End Proc: End CESS (DMT = Off)
Sequence: Seq01

Carol

Average
Step-wise
Yield : 98.8
Total bases = 35

A= 11, G= 4, C= 5, T= 15, 5= 0, 6= 0, 7= 0, 8= 0
(mixed bases= 0)

MW: 10706.0

5'> CCC ATT TAA ATA ATC TGA GCT TAT ATA TTT TGA GT <3'

Ethanol precipitate primers and set
up optimisation of PCR.

OPTIMISATION STRATEGY

	2mM	3mM	4mM	5mM
Buffer	35	35	35	35
MgCl ₂	14 μ l	21 μ l	28 μ l	35 μ l
dNTP	28 μ l	28 μ l	28 μ l	28 μ l
Primer	35 μ l	35 μ l	35 μ l	35 μ l
W-1	1.4	1.4	1.4	1.4
Taq	1.4	1.4	1.4	1.4
H ₂ O	207.2	207.2	207.2	207.2

2 Samples - 3 temperatures + 1 H₂O
= 7 tubes

$$A_{260} = 0.2411$$

$$\left(A_{260} \times \text{dil}^n \times 0.033 / \text{MW} \times 10^6 \right) = \text{Dir factor}$$

$$\left(\frac{0.2411 \times 200 \times 0.033}{10706} \right) \times 10^6 = \text{Dilution factor}$$

159 mg/ml


0.148 mm

1:7.43


RUNNING ORDER:

- 1) QX Hae III
- 2) 255 } 1
- 3) 256 } 4
- 4) 255 } 1
- 5) 256 } 4
- 6) 255 } 1
- 7) 256 } 4
- 8) 255 } 1
- 9) 256 } 4
- 10) 255 } 1
- 11) 256 } 4
- 12) 255 } 1
- 13) 256 } 4
- 14) H₂O - 1
- 15) - 2
- 16) - 3
- 17) - 4

Run 10ul on a
1% Agarose gel



- See p129 for
Running order



all Magnesium concentrations worked at all
 3 temperatures! Do 30 x 50 μ l PCR's
 on $\left. \begin{array}{l} 10 \times 1,1 \\ 10 \times 1,2 \\ 10 \times 2,2 \end{array} \right\}$ for Taq Polymorphism
 ie 31 tubes

Reagent	Volume	
Buffer	155 μ l	
MgCl ₂	62 μ l	
dNTP	124 μ l	4 μ l template
Primer	77.5 each	+ 46 μ l master mix
W-1	6.2 μ l	
Taq	6.2 μ l	Carry out reaction
H ₂ O	917.6	@ 56°C
	1426	

Digestions $\left. \begin{array}{l} 2.5 \\ 0.5 \end{array} \right\}$ μ l enzyme per tube
~~3660~~ tubes

$\left. \begin{array}{l} 7.5 \mu\text{l enzyme} \\ 60 \mu\text{l buffer} \end{array} \right\} 2.25 \mu\text{l}$
 $\left. \begin{array}{l} 15 \mu\text{l enzyme} \\ 120 \mu\text{l buffer} \end{array} \right\}$
 20 μ l

leave to digest overnight @ 37°C

Run samples on a 9% polyacrylamide gel
 @ 200V for 30 minutes

Samples used:-

Gel 1

112 - ^{TAG}1,1 - CC
 113 - 1,1 - -
 114 - 2,2 - GG
 115 - 1,1 - CC
 116 - 1,2 - GC
 117 - 1,1 - CC
 118 - 1,2 - GC
 119 - 1,1 - -
 120 - 1,1 - CC
 121 - 2,2 - GG
 122 - 1,2 - GC
 123 - 1,2 - GC
 124 - 1,2 - GC
 125 - 1,2 - GC
 126 - 2,2 - GG

Gel 2



127 - ^{TAG}1,1 - CC
 128 - 1,2 - GC
 129 - 1,1 - CC
 130 - 1,2 - GC
 132 - - CC
 133 - - GC
 134 - 1,2 - GG
 135 - 1,2 - GC
 136 - 1,1 - CC
 137 - 1,2 - GC
 138 - 1,1 - CC
 139 - 2,2 - GG
 140 - 1,1
 141 - 1,1
 142 -

12 x 1,1

11 x 1,2

4 x 2,2

Thus it appears that Allele (G) of my polymorphism is 100% associated with allele 2 of alisons Tag - and C is 100% linked to allele 1



FUTURE PLANS FOR tB845 POLYMORPHISM

- Fusions of AU rich region to β globin

Normal mRNA $\frac{1}{2}$ life
= 17hrs

- Band shift assays with both mRNAs
incubated with nuclear extract.

1st - ASTA using Adcel's Taq-man System


1/3/96

Set up Taq PCR again

~~10~~ x 100 μ l PCR'S > 4 MITS
Reactions

REAGENT	VOLUME
Buffer	40
MgCl ₂	40
dNTP	32
W-1	2
Taq	2
Primer	40
H ₂ O	244

- add Bul template
for each: -



- Obviously this has
not worked - why am
I getting such smears?

CONTROL EXPTS

- 1) Band extract DNA
before putting over
the column
- 2) Try Re-PCR'ing the
PCR product to
see if I get the
smears then
- 3) Run on acrylamide
& stain for protein

8845
1941
6904
20 each

Plan Expts

- PCR-up

1,1

2,2

1,2


} check that they are same for mice

- How much blood?

- Time scale?

IL1B peaks @ 6 hr

+ 18-24 hr



37° 15'

- Anti CD3 / CD19 beads - use these to pull off pure cell types

na of B cells

leave monocytes

Taqman 5 / 10 / 15 / 20 cycles

- Standard curve

Quantitation of steady state mRNA levels

⇒ Stability or accumulation of mRNA

eosinophils - acidic - stain eosin IL5

basophils - basic

neutrophils

mononuclear lymphocytes - from lymph - nucleus almost complete cell contains dense chromatin
monocytes - are large kidney-shaped nucleus much less dense than lymphocytes - because DNA is not dense

100 cells - 30 neutro
- 50 lympho
- mono

5/3/96


Do Taq PCR for further MUT-5 expts
- as previously set out (P115)



RUNNING ORDER

- 1) QX1174 Hae III
- 2) Taq
- 3) Taq
- 4) H₂O control

Also F5-BB₂ PCR for sequencing of 14123'



Running Order

- 1) ~~QX1174 Hae III~~
- 2) 165
- 3) 174
- 4) 149
- 5) 115
- 6) 157
- 7) 152
- 8) 139
- 9) 126
- 10) 121
- 11) 144
- 12) H₂O

DO 40 x 25 μ l reactions to confirm linkage of + ?? to Tag polymorphism

Taq PCR

3 samples from muts
1 positive control \Rightarrow 20 μ l reactions
1 H₂O control

REAGENT	VOLUME
Buffer	8
Mgcl ₂	8
dNTP	12.8
Primer	8 μ l
Taq	0.8
W-1	0.8
H ₂ O	97.8

RUNNING

- 1) 0.5 X Hae III
- 2) PVP 50
- 3) PVP 100
- 4) TWEEN
- 5) +ve control
- 6) H₂O control

Do ⁴⁸ 96 x 25 μ l PCR of new Polymorphism PCR

REAGENT	VOLUME
Buffer	120
MgCl ₂	96.48
dNTP	96
Primer	4.20
Taq	4.8
WT	4.8
H ₂ O	710.4
	110.4

2 μ l template per reaction

* Using the plate, some of the samples have evaporated.

- For digestion do 1st with ~~10~~ 10 μ l of sample digested with 0.5 μ l enzyme overnight and a further 1st with 0.25 μ l enzyme overnight.

mm1

6 μ l enzyme

18 μ l buffer

mm2

3 μ l enzyme

18 μ l buffer

7-19 0.5 μ l enzyme / 10 μ l

21-36 0.25 μ l enzyme / 10 μ l

6/3/96

Re-do FS-B62 PCR's on samples:-

152

139

126

121


144

REAGENT	VOLUME
Buffer	130
MgCl ₂	52
dNTP	104
Primer	< 65ul each
W-1	53 ul
Taq	52 ul
H ₂ O	769.6
TOTAL	1196

Add 8ul template
to each PCR.

RUNNING ORDER:-

- 1) QX HaeIII
- 2) 126
- 3) 139
- 4) 144
- 5) 152
- 6) 121
- 7) H₂O
- 8) 115 SS



7/3/96

D₀ 2 10 + 1440 50µl PCR_S for new polymorphism:-

REAGENT	VOLUME
Buffer	55 µl
MgCl ₂	22 µl
dNTP	44 µl
Primer	55 µl 27 µl each
W-1	2.2
Taq	2.2
H ₂ O	330 µl

4µl DNA per reacⁿ

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 95.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec


step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h49m21s



set up digests o/n.


5 - 0.5 μ l enzyme / 10 μ l


5 - 0.25 μ l enzyme / 10 μ l

Reaction mix 1: 5 μ l enzyme + 15 μ l buffer

2: 2.5 μ l enzyme + 15 μ l buffer

leave to cleave @ 37°C overnight \Rightarrow Run on
9% polyacrylamide gel.





Unable to cleave the top band
even with 10 μ l PCR product and
0.5 μ l enzyme @ 37°C overnight :-

all samples run today were
in equilibrium with +3953 polymorphism!

7 - 22 - GG

8 - 11 - CC

9 - 1,2 - GC

10 - 1,1 - CC

12 - 1,1 - CC

13 - 1,2 - GC

14 - 1,1 - CC


15 - 1,2 - GC

16 - 1,2 - GC

13/3

OPTIMISATION OF SCREENING FOR +6911

- dilute 7 μ l PCR product to 30 μ l with H_2O to cut out effects of buffers in PCR on enzyme.
- Add 3.5 μ l of buffer to each tube
- + 2 μ l enzyme. leave @ 37°C o/n.



Run 1/2 leave other 1/2 overnight for another night!



18/3/95



Do PCR's on 10 x F.B.

REAGENT	VOLUME
BUFFER	500
dNTP	400
MgCl ₂	200
Primer	250 each - 100µl biotinylated
Taq	20
w1	20
H ₂ O	2960
	4600µl


Redo Taq

~~P3/P2~~


- Running order.
- 1) ex haet
 - 2) 153
 - 3) 126
 - 4) 114
 - 5) 144
 - 6) 139
 - 7) 121



SQUIGGLES of: I11bb1t.connect March 18, 1988 17:07 J-R-U
 FOLDRNA of: I11bb1t Check: 877 from: 1 to: 871 March 7, 1988 17:34
 Length: 871 Energy: -228.5





G



MUTANT

SOURCE of: I11bb1t.connect March 18. 1998 20:57 S-5
 FOLDRNA of: I11bb1t Check: 708 from: 1 to: 971 March 18. 1998 20:12
 Length: 971 Energy: -228.8



22/3/96

Take Martin's IL-1 β probe in pUC 18 plasmid
and plate out on Amp/met + IPTG + XGal. Grow
@ 37°C overnight - Take one white colony
& inoculate 3 ml LB culture
Ampicillin
Meth.

LB RECIPE 1 LITRE

10g Bacto trypton
5g Yeast extract
5g NaCl

Autoclave before use.

3 ml culture was used to inoculate
a large 1 litre culture (2 x 500 ml). Again
ampicillin and methicillin were used.

Plasmid DNA is extracted by PEG
precipitation to give good, pure yields.
Using Ruth Herberts Protocol.

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☒ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.