


Express Mail Label No.: EV475172226US

Date of Deposit: February 9, 2005

Attorney Docket No. 24299-502



## IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICANTS: Duff *et al.*

APPLICATION NUMBER: 09/247,874

EXAMINER: Schnizer, Richard A.

FILING DATE: February 10, 1999

ART UNIT: 1635

FOR: THERAPEUTICS AND DIAGNOSTICS BASED ON A NOVEL IL-1B MUTATION

**MAIL STOP AF**  
Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

**DECLARATION OF DR. FRANCESCO S. DI GIOVINE UNDER 37 CFR 1.132**

Sir:

I, Francesco S. di Giovine, born on June 4, 1956, a University Senior Lecturer at the University of Sheffield, do hereby declare that:

1. I am one of the inventors of the above-identified application entitled "Therapeutics and Diagnostics Based on a Novel IL-1B Mutation."
2. I have an M.D. degree in Medicine and Surgery from the University of Florence (Italy), Faculty of Medicine, 1982, and a PhD degree in Molecular Immunology from the University of Edinburgh (UK), 1988, Faculty of Medicine.
3. I am actively engaged in researching genetic predispositions to various inflammatory diseases.
4. Working under my direction, members of my laboratory in early 1996 discovered and sequenced the IL-1B allele having a "C" rather than a "G" at the position corresponding to +6912 in Figure 1, which is nucleotide 8845 of SEQ ID NO: 1. This allele is also referred to as IL-1B (+6912) allele 2.
5. A copy of the laboratory notebook ("the Notebook") kept by my assistant Carol Campbell and reviewed by me on an ongoing basis is attached hereto as Appendix A. This notebook contains experiments performed in my laboratory from August 24, 1995 to March 22, 1996. The page numbering is located in the lower left or right-hand corner of each sheet. I confirm that this is a true and complete copy of the notebook kept by Ms. Campbell during the time-period in question.

Applicants: Duff *et al.*  
U.S.S.N.: 09/247,874

6. As described in the instant application (Example 1, pp. 36-37), a PCR product corresponding to the 3'UTR of the IL-1B gene was amplified from human genomic DNA and sequenced. The primers used in the amplification were disclosed in the application as SEQ ID NO: 3 (5'-GTCCCCACATTCTGATGAGCAAC-3') and SEQ ID NO: 4 (5'-TGCAGCACTCAGCAATGAGGAG-3'), which bind to regions of the IL-1B gene corresponding to positions +6720 to +6742 and +7102 to +7123, when the IL-1B gene is numbered in accordance with the numbering of Figure 1. These primers were designed on August 29, 1995 and are represented as primers F<sub>2</sub> and B<sub>1</sub> on page 3 of the Notebook. The determination that a G to C change at the +6912 location was performed by sequencing of the fragment amplified by the F<sub>2</sub> and B<sub>1</sub> PCR primers, the result of which was recorded on February 21, 1996 on page 116 of the Notebook. The oligonucleotide primer used for sequencing of the +6912 allele 2 is shown on page 117. A sequence-specific oligonucleotide primer hybridizing to nucleotides +6913 to +6947 of the IL-1B gene is described on page 118. The location of the G to C change is further evidenced on page 124, entitled "Further plans for the +8845 polymorphism," because as stated above, nucleotide 8845 of SEQ ID NO: 1 corresponds to the IL-1B (+6912) location.

7. I understand that the Examiner, while admitting that we have discovered a G to C polymorphism at position +6912, stating that neither the specification or my previously filed declaration supports a 9721 nucleotide sequence of IL-1B with a C at position +6912. (See, Office action at page 5). The Examiner states that "[s]ince polymorphisms can occur throughout a molecule, one cannot assume that there are no other polymorphisms linked to position +6912 within the 9721 bases of the IL-1B gene, and that the sequence of the rest of the 9721 nucleotides is identical to that reported in the prior art." (See, Office action at pages 5-6).

8. I believe that the C polymorphism at position +6912 is a single nucleotide polymorphism, or "SNP." A SNP is known in the art as a DNA sequence variation among individuals in which the purine or pyrimidine base (as guanine) of a single nucleotide in the genome has been replaced by another such base (as cytosine). Therefore, I believe that one of skill in the art would recognize that the sequence of the rest of the 9721 nucleotides of IL-1B is identical to the wild-type sequence, which has been reported in the art, and that upon identifying the SNP at position +6912, it was unnecessary to re-sequence the entire IL-1B gene.

9. Regarding the Examiner's statement regarding linkage analysis of the IL-1B +6912 with other IL-1B polymorphisms, I note that genetic linkage between alleles of a given gene such as IL-1B does not indicate the presence of multiple sequence variations in the IL-1B gene of an individual. Rather, I believe that linkage analysis is useful to demonstrate co-segregating polymorphisms that contribute to a given disease or disorder. As shown in the Notebook, under my direction, linkage analysis was performed on or about February 23, 1996 between the IL-1B (+6912) locus and the IL-1B *taq* locus, and the result is shown on page 129 of the Notebook. This analysis demonstrates that the IL-1B (+6912) allele 1 (termed "Allele (G)" in the notebook) is 100% associated with the IL-1B *taq* allele 2, while the IL-1B (+6912) allele 2 (termed "Allele (C)" in the notebook) is 100% associated with the IL-1B *taq* allele 1, further

Applicants: Duff *et al.*  
U.S.S.N.: 09/247,874

demonstrating that my laboratory had identified a novel allele having a "C" at position +6912 of the IL-1B gene.

10. At the time of these experiments, the human IL-1B gene sequence published by Clark *et al.* (Nucleic Acids Research 14(20):7897-7914 (1986)) was regarded as the standard sequence for human IL-1B. This sequence is also deposited in GenBank under the accession number X04500. The Clark *et al.* sequence shows a "G" at position +6912, in contrast to the allele we discovered, which contains a "C" at position +6912. We named the "G" variation "allele 1" and the "C" variation "allele 2." The statement in the application indicating that the IL-1B allele 1 has a cytosine at position +6912 and that allele 2 has a guanine at that position is a typographical error. I believe that one of skill in the art would recognize the existence of this error based on the teachings of the application and, further, in view of the contents of the Notebook. Our measurements of allele frequency presented in the patent application (e.g., Example 2, pp. 37-38) demonstrate that allele 1 is the more frequent allele and may therefore be considered the wild-type allele. Moreover, multiple publications and database entries have presented the Clark *et al.* nucleic acid sequence as the wild-type sequence. (See, e.g., US Patent numbers 5,686,246; 6,720,141; 6,730,476; and 6,746,839; and GenBank Accession number P01584).

11. In conclusion, in February, 1996, prior to the time of filing of the present application, I was in possession of a novel IL-1B sequence with a "C" at position +6912 from the transcription start site, which is position 8845 of a nucleic acid sequence numbered in accordance with SEQ ID NO: 1. This allele is termed the IL-1B (+6912) allele 2 and is substantially less common than the IL-1B (+6912) allele 1 known in the prior art.

Applicants: Duff *et al.*  
U.S.S.N.: 09/247,874

12. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

*Francesco di Giovine*

---

FRANCESCO S. DI GIOVINE, M.D., Ph.D.

Dated: January 31<sup>st</sup>, 2005

TRA 1988391v1

Graham

2855-322

Perkin elmer. seq

User No  
13

## PCR lab.

23

Mohammed - Work 0116 254 1414 bleep  
4099  
- Home 0116 232 2675

franco Hotel Fax:- 03 9622 8877

10 Mrs diff - Fax @

12pm  $\Rightarrow$  10am

Niceday Manuscript Books

## Niceday Manuscript Books

## Range of Sizes and Formats

Size	Format	No. of Sheets	Product Code
<b>Everyday Range</b>			
A6	Feint & indexed A-Z	96	No.237204
A5	Feint & indexed A-Z	96	No.217539
A4	Feint & indexed A-Z	96	No.215221
A6	Feint Ruled	96	No.118080
A5	Feint Ruled	96	No.144300
A4	Feint Ruled	96	No.142610
A6	Feint & indexed A-Z Wirebound	80	No.287533
A5	Feint & indexed A-Z Wirebound	80	No.178695
A4	Feint & indexed A-Z Wirebound	80	No.164569
A6	Feint Ruled Wirebound	80	No.199230
A5	Feint Ruled Wirebound	80	No.179273
A4	Feint Ruled Wirebound	80	No.233139

### **Environment Range (100% Recycled)**

A6 Feint Ruled 96 No.174951  
A5 Feint Ruled 96 No.283590  
A4 Feint Ruled 96 No.118367

### **Economy Range**

A6	Feint Ruled	96	No.135969
A5	Feint Ruled	96	No.174741
A4	Feint Ruled	96	No.240288

GRANT NO. [56833]

(m) PNAS - 1993 ao 6 2295.

(D) J. Immun 94 153-2 712.

(m) ADVANCES IN NEUROIMMUNOLOGY 92 2 No 1.




D) Int J. immuno. pharm. 92 16 No 3

m) Mol & Cel Biol 95 15 1 112-119.

(M) Mol Immun 1995 32 No 8 541-


24/8/95

caat 46 bp repeats containing 36' recognition site



+3953 (C-T)

5 x 46 bp VNTR



-889 (C-T)

Primer 1: IL1B - Annealing temperature 55.9°

Primer 2: IL1A - Annealing temperature 53.4°

IL1B

PmALPS stimulation - 2982 - 2795

LPS stimulation - 3757 - 2729

NFKB sites (within promoter) - 296/-286  
- 2761/-2753

29/8/95

## Macrophage migration Inhibitory Factor (mIF)

- Exists preformed in macrophages & monocytes
- mIF inhibits TNF secretion
- Inhibits stimulated NO secretion by macrophages.
- 90% identity between mouse and human.
- mIF secretion is induced by TNF $\alpha$ , IFN Gamma
- $\leq 1\text{Kb}$  - 3 exons separated by 2 introns of 189 & 95bp
- Multiple genes in mouse but not human.

### Primers - IL-1 $\beta$ 3'

Primers were designed using macvector. These primers span 3' AU rich region of IL-1 $\beta$ . The AU rich region confers instability upon an mRNA sequence - thus cytokines do ~~not~~ have them in order to return cytokine levels to norm as soon as an attack is over (TATTAT.)

F<sub>1</sub> 5' CAAG CAGAAAA CATGCC GTC 3' (Tm 57.0)

F<sub>2</sub> 5' GCTCCC ACA TTCT GATG AGCAAC 3' (Tm 57.4)

F<sub>3</sub> 5' CAATT GATGAG CAA CGGCC 3' (Tm 56.2)

B<sub>1</sub> 5' TGCAAGCACT CAGCAA TGAGGAG 3' (Tm 57.8)

30/8/95

Primers spanning the A-T region of IL1 $\alpha$  were also designed

F4 -  $^5'$  ATA GCA TAA GTT TCT TGG ACC TCA G.  $^3'$

B3/ccc/IL1 $\beta$  -  $^5'$  CAG ATA CTG GAAAA CC A GGC GTAGG  $^3'$

B2/ccc/IL1 $\alpha$  -  $^5'$  GCT TGT AGG ACT TGA TTG CAGGTG C  $^3'$

F4/IL1 $\alpha$ /ccc + B2/IL1 $\alpha$ /ccc give a 522 bp fragment.

F5/IL1 $\alpha$ /ccc + B2/IL1 $\alpha$ /ccc give a 984 bp fragment.

Primers were designed using 'macvector' and checked using BLAST searching to ensure that they matched no other sequence in the human genome.

**GIBCO BRL**

**Taq DNA Polymerase**

**Cat. No. 18038-026**

**Lot No. FET401** 500 units; 5 U/ $\mu$ l

Exp. Date: 05/97. Store at -20°C (not frost-free).

**LICENSED FOR PCR**

**Description:**  
Taq DNA Polymerase is isolated from *Thermus aquaticus* YT1. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kD. Taq DNA polymerase is heat-stable and will synthesize DNA at elevated temperatures from single-stranded templates in the presence of a primer. Effective January 1, 1994, Life Technologies™ has modified the unit assay for Taq DNA polymerase, effectively doubling the enzyme concentration for consistent PCR performance.

**Components:**

18038-026	Taq DNA Polymerase	Lot No. <b>FET401</b>
Y02028	10X PCR Buffer	Lot No. <b>FDJ102</b>
Y02016	50 mM Magnesium Chloride	Lot No. <b>FCC102</b>
90238	1% W-1	Lot No. <b>FET403</b>

**Unit Definition:**  
One unit incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

**Storage Buffer:**  
20 mM Tris-HCl (pH 8.0)  
0.1 mM EDTA  
1 mM DTT  
50% (v/v) glycerol  
Stabilizers

**Unit Assay Conditions:**  
25 mM TAPS (pH 9.3)  
50 mM KCl  
2 mM MgCl<sub>2</sub>  
1 mM DTT  
0.5 mg/ml activated salmon sperm DNA  
0.2 mM dATP, dCTP, dGTP, dTTP

**10X PCR Buffer:**  
200 mM Tris-HCl (pH 8.4)  
500 mM KCl

**The PCR Buffer supplied as a 10X concentrate should be diluted for use.**

The 1% solution of the detergent W-1 can be added at a final concentration of 0.05% (v/v) and may improve the thermostability of the enzyme. Store solution at -20°C and thaw at 37°C before use.

Doc. Rev.: 022895

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE IN HUMANS OR ANIMALS.

Printed on Recycled Paper

LIFE TECHNOLOGIES

## P.C.R.

dNTP mix.      10 $\mu$ l dATP  
                  10 $\mu$ l dCTP  
                  10 $\mu$ l dGTP  
                  10 $\mu$ l dTTP  
                  360 $\mu$ l H<sub>2</sub>O       $\Rightarrow$  10mM working stock  
                  400 $\mu$ l

P.C.R. buffer - supplied as 10x.  $\Rightarrow$  working conc<sup>n</sup> = 1x.

1 $\mu$ l detergent per 1 $\mu$ l Tag

INITIAL DENATURE - 96° - 2'  
DENATURE 94° - 1'      }  
56°C Anneal 1'      } X 35.  
ELONG 72° 1'      }

FINAL ELONG 72° 5'  
4°

D/S DNA	1 O.D. = 50 $\mu$ g/ml
RNA	1 O.D. = 40 $\mu$ g/ml
S/S DNA	1 O.D. = 33 $\mu$ g/ml

P.C.R. 1L-1B 3' END

<u>Reagent</u>	<u>Stock</u>	<u>use</u>	<u>final</u>
10 x P.C.R. buffer	-	5μl	1 x
MgCl <sub>2</sub>	50mM	2.5μl	2.5mM
dNTP's mix	10mM (each)	4μl	0.2mM
Taq Pol (GIBCO)	50/μl	0.2μl	10
W-1 (detergent)	-	0.2μl	-
Template	50μg/ml	2μl	100ng/reac
Primer mix	2μM each	5μl	1μM
H <sub>2</sub> O	-	31.1	-
		50μl	
2 drops mineral oil			

F<sub>1</sub> + B<sub>1</sub> annealing Temperature = 56.0°

F<sub>2</sub> + B<sub>1</sub> annealing Temperature = 53.3°

P.C.R. CYCLES

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec } X 35

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h53m59s

If a P.C.R. Reaction does not work first time - alter  
MgCl<sub>2</sub> concn or annealing Temperature or

1/9/95

Primers F1/IL-1B/ccc  
F2/IL1B/ccc  
B1/IL1B/ccc

were made and supplied  
in dry form. They were  
extracted and purified  
by ethanol pptn.

MM 737 = B1/IL1B/ccc

MM 738 = F2/IL1B/ccc

MM 739 = F1/IL1B/ccc

	B1	F1	F2
1) Resuspend oligo in 200 $\mu$ l H <sub>2</sub> O	✓	✓	✓
2) Place in eppendorf	✓	✓	✓
3) Add 10 $\mu$ l 3M NaAc.	✓	✓	✓
4) Add 300 $\mu$ l 69% EtOH	✓	✓	✓
5) -70°C 1hr. - 11:30am	✓	✓	✓
6) Spin 15' 12,000 rpm	✓	✓	✓
7) Remove supernatant	✓	✓	✓
8) Wash in 70% ETOH	✓	✓	✓
9) Remove supernatant	✓	✓	✓
10) Dry in vacuum ~10'	✓	✓	✓

\* When you release vacuum - turn black  
mark away from hole - ie: off.

Remove hose & turn off H<sub>2</sub>O first

\* Otherwise, you will flood your samples

\* AFTER DRYING THE OLIGO SHOULD BE RESUSPENDED IN 50 $\mu$ l \*

O.D

Programme 10 - Walborg.

Step - 10.

press RS

Recall press O ⇒

insert blank

H<sub>2</sub>O  
flush

7

Put tube in touch fill press RS

O.D. CALCULATION-

add 170  $\Rightarrow$  200  $\mu\text{l}$ . - take 5  $\mu\text{l}$  + 995  $\text{H}_2\text{O}$

$A_{320}$  low. -

$$\underline{B_1} \quad A_{260} \quad 0.3188 \quad = \quad 3$$

$$= [J] = 10.52 \text{ mg/ml}$$

$$\text{MW} = 6787.4$$

$$6787.4 \text{ g/l} \text{ litre} = 1 \text{ M}$$

$$6784 \text{ mg/ml} = 1 \text{ M}$$

$$10.52 \text{ mg/ml} = 1.55 \text{ mM}$$

$$= 1.55 \mu\text{M} \quad \Rightarrow \text{ dilute 1:77.5}$$

$$1 \mu\text{l} + 77 \mu\text{l} \text{ H}_2\text{O}$$

$$400 \mu\text{M soln.} \quad \Rightarrow \quad 1 \text{ ml} + \text{ soln.} = 195 \mu\text{l} \text{ H}_2\text{O} + 65 \mu\text{l} \text{ primer}$$

Primer mix	= 4 $\mu\text{l}$ $B_1$ + 4 $\mu\text{l}$ $F_1$	= 260 $\mu\text{l}$ stock
	+ <del>76</del> 152 $\mu\text{l}$ $\text{H}_2\text{O}$	

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
$B_1$	1.0000	0.0191	0.2098	0.3188	0.6362	1.5718	11.988
$F_1$	2.0000	-0.004	0.1123	0.1435	0.7863	1.2685	5.0951
$P_2$	3.0000	-0.025	0.1237	0.1154	1.0593	0.9440	3.4762
$\text{H}_2\text{O}$	4.0000	-0.020	0.0000	0.0934	0.1770	5.6494	6.4125

Column 1

14:13:08 ,30/ 8/95

Run ID : MM 737

Cycle : 002 UMOL  
End Proc: End CESS (DMT = Off)  
Sequence: Seq01

vector (1000)

Average

Step-wise

Yield : 98.4

Total bases = 22

A= 7, G= 7, C= 5, T= 3, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 6787.4

5'> TGC AGC ACT CAG CAA TGA GGA G <3' B1 /IL1B /ccc  
a056

a055

Column 2

T= 3

14:13:09 ,30/ 8/95

Run ID : MM 738 C= 5  
Cycle : 002 UMOL  
End Proc: End CESS (DMT = Off)  
Sequence: Seq02

40  $\mu$ M

$\Rightarrow$  Stock Soln. (400  $\mu$ M)

Average

Step-wise

Yield : 96.8

Total bases = 23

A= 6, G= 4, C= 8, T= 5, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 6965.6

5'> GCT CCC ACA TTC TGA TGA GCA AC <3' F2 /IL1B /ccc

Column 3

14:13:09 ,30/ 8/95

Run ID : MM 739

Cycle : 002 UMOL  
End Proc: End CESS (DMT = Off)  
Sequence: Seq03

1

546  $\mu$ M

1:4 soln.

91  $\mu$ l H<sub>2</sub>O = 156  $\mu$ l soln.

Average  
Step-wise  
Yield : 96.8  
Total bases = 21  
A= 8, G= 4, C= 7, T= 2, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 6390.2

5'> CAA GCA GAA AAC ATG CCC GTC <3' F1 /IL1B /ccc

1:20 diln.

H<sub>2</sub>O

F

$$A_{260} = 0.1435 \times \text{dil factor (1000)}$$

$$= 143.5 \times 33$$

$$= \text{OD of } \underline{4.735}$$

1446

3550.4

2336.

3438.4

10770.8

$$M_w = 6390$$

$$6390 \text{ g/l} = 1M$$

$$6390 \text{ mg/ml} = 1M$$

$$\cancel{6390 \text{ mg/l}} = 1M$$

$$4.735 \text{ mg/ml} = 700 \mu M$$

$$= 65 \mu \text{l primer} + 130 \mu \text{l H}_2\text{O} \quad \xrightarrow{1:2} \text{Stock Soln. (400} \mu \text{M})$$

F2

$$A_{260} = 0.1154 \times 1000 =$$

$$115.4 \times 33 = 3808.2$$

$$= 3.808 \text{ mg/ml}$$

$$M_w = 6965$$

$$6965 \text{ g/l} = 1M$$

$$6965 \text{ mg/ml} = 1M$$

$$3.808 \text{ mg/ml} = 546 \mu M$$

$\Rightarrow$  to get  $400 \mu M$  soln  $1:1.4$  soln

~~65~~  $65 \mu \text{l primer} + 91 \mu \text{l H}_2\text{O} = 156 \mu \text{l soln}$

To get from  $400 \rightarrow 20 \mu M$   $1:20$  dil.

Primo mix =  $2 \mu \text{l}$  <sup>each primer</sup> +  $76 \mu \text{l}$  of  $\text{H}_2\text{O}$

$$P_1 = B_1 + F_1 \quad \times 4$$

$$P_2 = B_1 + F_2$$

## IL-1B TEST PCR

Two PCR reactions were set up 1 with F1/IL1B/ccct + B1/IL1B/ccc (803 bp fragment) and 1 with F2/IL1B/ccc + B1/IL-1B/ccc (403 bp fragment)

Template used = 3 x normal genomic samples

⇒ Set up 8 reactions

Master mix:- Buffer = 40  $\mu$ l ✓

MgCl<sub>2</sub> = 20  $\mu$ l ✓

dNTP's = 32  $\mu$ l ✓

Taq Pol = 1.6  $\mu$ l ✓

W1 = 1.6  $\mu$ l ✓

Template

$\ddagger$  H<sub>2</sub>O 24.8  $\mu$ l ✓

344.0  $\mu$ l

Pipette 43  $\mu$ l of Master mix into each tube ✓

into tubes 1, 2, 3 and 4 add 5  $\mu$ l Primer mix 1 ✓

into tubes 5, 6, 7 and 8 add 5  $\mu$ l Primer mix 2

tubes 1 and 8 are H<sub>2</sub>O controls

"lib 02:1" May 25 2004 now 2:00 p.m.

0.5  $\mu$ l 20  $\mu$ l 5  $\mu$ l + 10  $\mu$ l = 1.0  $\mu$ l dilution

1.0  $\mu$ l + 1.0  $\mu$ l = 2.0  $\mu$ l

1.0  $\mu$ l + 1.0  $\mu$ l = 2.0  $\mu$ l

TUBE NO.\*DNA USED

German controls

1/	- H <sub>2</sub> O control	+ Primer mix 1	H <sub>2</sub> O
2/		+ Primer mix 1	Sample No 1
3/		" "	" " 3
4/		" "	" " 7
5/		+ Primer mix 2 (oops Primer 1 added by mistake)	" " 8
6/		" "	1
7/		" "	3
8/	- H <sub>2</sub> O control	+ Primer mix 2	H <sub>2</sub> O

Tubes are 'loaded' as above - before adding DNA -  
Wax pellets are added to each tube they are  
then melted @ 75°C for 5' and cooled to 4°C  
for 3' DNA samples (2  $\mu$ l) are then added  
on top of the wax layer before putting on  
the thermocycler.

No 61 - 75° 5' } wax melting  
62 4° 3' }

\* The volume of reagent which should be put on top of  
the wax layer should be 10-30  $\mu$ l - To ensure equal  
mixing this can be achieved by adding Master mix to  
a DNA / H<sub>2</sub>O solution which is wax sealed &

IL-1 $\alpha$  PrimersDATA

- 394 Synthesis Setup Listing - (Version 2.00)

Column 1

13:50:56 , 31/ 8/95

Run ID : MM 740

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq01

GCT TGT AGG ACT TGA TGC

AGG TGG 3'

Average

Step-wise

Yield : 98.5

Total bases = 25

A= 4, G= 10, C= 3, T= 8, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 7775.0

WLS

11420

5'> GCT TGT AGG ACT TGA TTG CAG GTG G <3' B2/ccc/IL-1 $\alpha$ .

ATA GCA TAA GTT TCT GGG ACC TCA

ATA CTG GAA AAC CAG GCG TAG

Column 2

13:50:57 , 31/ 8/95

Run ID : MM 741

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq02

Average

Step-wise

Yield : 98.8

Total bases = 25

A= 7, G= 6, C= 5, T= 7, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 7676.0

11420

5'> ATA GCA TAA GTT TCT GGG ACC TCA G <3' F4/ccc/IL-1 $\alpha$ .

52.7°  $\Rightarrow$  522 bp fragment

54.6°  $\Rightarrow$  984 bp fragment

Column 3

13:50:57 , 31/ 8/95

Run ID : MM 742

Cycle : 002 UMOL

End Proc: End CESS (DMT = Off)

Sequence: Seq03

Average

Step-wise

Yield : 98.2

Total bases = 25

A= 9, G= 8, C= 5, T= 3, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 7742.0

5'> CAG ATA CTG GAA AAC CAG GCG TAG G <3'

F5/ccc/IL-1 $\beta$

IL-1 $\alpha$  Primers

DATA

- 394 Auto Analysis Listing -

Time: 16:18:39 , 1/ 9/95

Column £1	Column £2	Column £3			
Seq: Seq01 B2	Seq: Seq02 F4	Seq: Seq03 F5			
Overall: 69.9	Overall: 75.1	Overall: 66.0			
ASWY: 38.5	ASWY: 98.8	ASWY: 98.3			
Num Base	ASWY	Num Base			
		ASWY			
2 G	100.0	2 A	100.0	2 G	100.0
3 T	97.1	3 C	92.5	3 A	97.3
4 G	94.9	4 T	94.9	4 T	97.7
5 G	96.2	5 C	96.2	5 G	97.6
6 A	96.9	6 C	96.9	6 C	96.2
7 C	95.9	7 A	97.4	7 G	96.8
8 G	96.4	8 G	97.8	8 G	96.7
9 T	96.9	9 G	98.1	9 A	96.6
10 T	97.2	10 G	98.3	10 C	96.9
11 A	96.8	11 T	98.5	11 C	97.2
12 G	97.1	12 C	98.6	12 A	97.5
13 T	97.4	13 T	98.7	13 A	97.7
14 T	97.6	14 T	98.8	14 A	96.8
15 C	97.6	15 T	98.9	15 A	97.1
16 A	97.8	16 G	99.0	16 G	97.3
17 G	97.9	17 A	99.0	17 G	97.4
18 G	98.1	18 A	99.1	18 T	97.6
19 A	98.0	19 T	98.8	19 C	97.7
20 T	98.1	20 A	98.9	20 A	97.8
21 G	98.2	21 C	98.6	21 T	97.9
22 T	98.3	22 G	98.6	22 A	98.0
23 T	98.4	23 A	98.7	23 G	98.1
24 C	98.5	24 T	98.8	24 A	98.2
25 G	98.5	25 A	98.8	25 C	98.3

B2|ccc|IL-1 $\alpha$  = 5' G C T T G T A G G A C T T G A T G G A G G T G G 3'

F4|ccc|IL-1 $\alpha$  = 5' A T A G C A T A A G T T T C T G G G A C C T G 3'

F5|ccc|IL-1 $\alpha$  = 5' C A G A T A C T G G A A R A C C A G G C G T A G 3'

F4 + B2 - Annealing Temperature = 52.7°  $\Rightarrow$  522 bp fragment  
F5 + B2 - Annealing Temperature = 54.6°  $\Rightarrow$  984 bp fragment

1145

129  
123

4/9/95

Purify, by ethanol precipitation newly made primers  
 B2/L12/ccc, F4/L12/ccc, F5/L12/ccc

- 1) Resuspend oligos in 200 $\mu$ l Pure H<sub>2</sub>O.
- 2) Put 100 $\mu$ l into ~~one~~ fresh eppendorf.
- 3) Add 10 $\mu$ l 3M NaAc.
- 4) Add 300 $\mu$ l EtOH (100%).
- 5) Freeze @ -70°C ~1hr.
- 6) Spin, 12,000 rpm - 15
- 7) Remove SN
- 8) Wash: pellet in 70% EtOH
- 9) Spin 5' 12,000 rpm
- 10) Dry Under vacuum ~10'.
- 11) Resuspend 500 $\mu$ l H<sub>2</sub>O.
- 12) Measure O.D. (Prog 10. Warburg)

\* To put vacuum on - turn black mark away from 'hole'  
 put on tube and turn on top - to turn off, pull tube  
 off - then switch off H<sub>2</sub>O then lift nozzle \*

$$\begin{aligned}
 \text{B2} \quad A_{260} - 0.5413 \times \text{dil factor.}^{(200)} \\
 &= 108.26 \times 33 \\
 &= 00. = 3572.58
 \end{aligned}$$

$$\begin{aligned}
 \text{MW} = 7775 &= 7775 \text{ g/l} = 1 \text{ m} \\
 7775 \mu\text{g/l} \text{ u} &= 1 \text{ m} \\
 3572.58 &= 460 \mu\text{m}
 \end{aligned}$$

Dilution factor = 27

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
1.0000	-0.010	0.0104	0.0161	0.7798	1.2824	11.850	0.9111
32 2.0000	0.0630	0.3096	0.5413	0.5156	1.9397	20.491	21.209
F4 3.0000	-0.023	0.1182	0.2347	0.5475	1.8264	23.812	11.122
F5 4.0000	0.0264	0.1240	0.2105	0.5304	1.8853	12.133	8.0636

Run P.C.R. Samples on 1% Agarose Gel. Add 1 $\mu$ l EtBr to gel  
(P1011)  
Make small 30ml gel. 0.5g Agarose  
50ml TBE (1x)

- Take 10 $\mu$ l PCR product
- Add 4 $\mu$ l loading dye

also run 0.5-1 $\mu$ l  $\alpha$ X HaeIII marker

loading order.

GEL 1

- 1)  $\alpha$ X HaeIII ~~too~~ Molecular weight marker
- 2)  $\text{H}_2\text{O control}$  Sample No 1 (Primers B<sub>1</sub>/F<sub>1</sub>)
- 3) " 3
- 4) " 7
- 5) " 8
- 6)

Gel 2

- 1)  $\alpha$ X HaeIII Molecular WT marker
- 2) ~~H2O~~ Sample No 1 (Primers B<sub>1</sub>/F<sub>2</sub>)
- 3) No 3
- 4) H<sub>2</sub>O control

Gel was run @ 80V for ~ 30'

RESULT No bands were seen except those of primers - so re-plan expt. optimise conditions such as magnesium concentration and annealing Temperature. Increasing magnesium concentration lowers the specificity of the reaction as does lowering the annealing temp.

LIB PCR - MARK II

OPTIMISATION

alter -  $MgCl_2$

Annealing Temp.

$MgCl_2$  - was 2.5 use: 2.0 3.0 3.5

Annealing Temperature was 53° - try 50 48  
56° 53 50

Block 1 Annealing temp: 56°C

Samples 1, 3

Magnesium 2.0 3.0 3.5  $\mu$ l

Master mix :- 1 x 9

Buffer 5  $\mu$ l

$MgCl_2$  2.5  $\mu$ l (2.5 mM)

dNTP's 4  $\mu$ l

Taq 0.2  $\mu$ l

W1 0.2  $\mu$ l ✓

Template 2  $\mu$ l

Primer 5  $\mu$ l ✓

H<sub>2</sub>O 31.5

50

Mix 2

5  $\mu$ l

(3.5 mM)

3.5  $\mu$ l

4  $\mu$ l

0.2  $\mu$ l

0.2  $\mu$ l

2  $\mu$ l

5  $\mu$ l

29.5

50

Mix 3

5  $\mu$ l ✓ (4.5 mM)

4.5  $\mu$ l

4  $\mu$ l

0.2  $\mu$ l

0.2  $\mu$ l

2  $\mu$ l

5  $\mu$ l

29.5

50

31.5  
+ a  
✓  
28.5

30.5  
x 9  
27.5

29.5  
x 9  
26.5

	(1) mm1	(2) mm2	(3) mm3
① Block 1: 56°	$\frac{(10)(100)}{1, 2, \text{H}_2\text{O}}$	1, 2, $\text{H}_2\text{O}$	1, 2, $\text{H}_2\text{O}$
② Block 2: 53°	1, 2, $\text{H}_2\text{O}$	1, 2, $\text{H}_2\text{O}$	1, 2, $\text{H}_2\text{O}$
③ Block 3: 50°	1, 2, $\text{H}_2\text{O}$	1, 2, $\text{H}_2\text{O}$	1, 2, $\text{H}_2\text{O}$

### O.D Primers

$$F_4, A_{260} = 0.2347$$

~~$$0.2347 \times 200 \times 33 = \text{OD}$$~~

~~$$\text{O.D} = 1549.02$$~~

Dilution factor = 10

~~$$Mw = 76760 = 1\text{M} = 76760 \mu\text{g}/\mu\text{l}$$~~

~~$$\Rightarrow 1549 = 202 \mu\text{M}$$~~

~~$$F_5, A_{260} = 0.2105 \times 200 \times 33 = \text{OD}$$~~

~~$$\text{O.D} = 1389.3$$~~

$$Mw = 6390 = 220 \text{ mM}$$

F5 - dilution factor = 11  
to give 20 mM soln

$$\frac{(A_{260} \times \text{diln} \times 0.033 / \text{mw}) \times 10^6}{20} = \text{diln factor}$$

$$\frac{F_5 \times (1.3893) \times 10^6}{20} = 10$$

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h48m11s

Samples A

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h52m59s

Samples B

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 54.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h59m23s

Samples C

IL-1 $\alpha$  PCR

PROVISIONAL.  
IL 1 $\alpha$  PCR PROTOCOL



REAGENT	STOCK	USE	FINAL
10 X BUFFER		5 $\mu$ l	1X
MgCl <sub>2</sub>	50mM	2.5 $\mu$ l	2.5mM
dNTP's	10mM	4 $\mu$ l	0.2 mM
Taq Pol	5U/ $\mu$ l	0.2 $\mu$ l	1U
W-1		0.2 $\mu$ l	
TEMPLATE	50 $\mu$ g/ml	2 $\mu$ l	100ng/react
primer mix	20 $\mu$ M each	5 $\mu$ l	1 $\mu$ M
H <sub>2</sub> O		<u>31.1<math>\mu</math>l</u>	
		<u>50<math>\mu</math>l</u>	

CYCLES:

INITIAL DENATURATION	96°C	2MINS
DENATURATION	94 °C	1 MIN
ANNEALING	F4 (52 °) F5(54 °)	1 MIN
ELONGATION	72°C	X35 1 MIN
ELONGATION	72°C 4°C	5 MINS INFINITY

Primer mix 3 = F4/IL1 $\alpha$ /ccc + B2/IL1 $\alpha$ /ccc @ 20 $\mu$ M

Primer mix 4 = F5/IL1 $\alpha$ /ccc + B2/IL1 $\alpha$ /ccc @ 20 $\mu$ M



5/9/95

Run PCR samples on a 1% ~~Agarose~~ <sup>Agarose</sup> Gel

10 μl PCR products

1 μl bromophenol blue

load 14 μl

Loading order

1)	1 μl Q <sub>X</sub> Hae III Mw marker	18)	1 μl Q <sub>X</sub> Hae III
2)	14 μl A <sub>1</sub> (1)	19)	14 μl B <sub>3</sub> (1)
3)	" A <sub>1</sub> (2)	20)	" B <sub>3</sub> (2)
4)	A <sub>1</sub> H <sub>2</sub> O	21)	B <sub>3</sub> (H <sub>2</sub> O)
5)	A <sub>2</sub> (1)	22)	" C <sub>1</sub> (1)
6)	A <sub>2</sub> (2)	23)	" C <sub>1</sub> (2)
7)	A <sub>2</sub> (H <sub>2</sub> O)	24)	" C <sub>1</sub> (3)
8)	A <sub>3</sub> (1)	25)	" C <sub>2</sub> (1)
9)	A <sub>3</sub> (2)	26)	" C <sub>2</sub> (2)
10)	A <sub>3</sub> (H <sub>2</sub> O)	27)	" C <sub>2</sub> (H <sub>2</sub> O)
11)	B <sub>1</sub> (1)	28)	" C <sub>3</sub> (1)
12)	B <sub>1</sub> (2)	29)	" C <sub>3</sub> (2)
13)	B <sub>1</sub> (H <sub>2</sub> O)	30)	" C <sub>3</sub> (H <sub>2</sub> O)
14)	B <sub>2</sub> (1)	31)	" Q <sub>X</sub> Hae III
15)	B <sub>2</sub> (2)		
16)	B <sub>2</sub> (H <sub>2</sub> O)		
17)	1 μl Q <sub>X</sub> Hae III Mw marker		

Run gel at 120 v. for ~ 1/2 hour.

Set up 1L12 PCR as per provisional protocol P19.

MM(1) buffer = 15 $\mu$ l  
MgCl<sub>2</sub> = 7.5 $\mu$ l  
dNTP's = 2.12 $\mu$ l  
Taq Pol = 0.6 $\mu$ l  
W-1 (F1/F2) = 0.6 $\mu$ l  
Primer mix = 2.6 $\mu$ l  
H<sub>2</sub>O = 93.3 $\mu$ l  
150 $\mu$ l

MM(2)  
buffer 15 $\mu$ l  
MgCl<sub>2</sub> 7.5 $\mu$ l  
dNTP's 12 $\mu$ l  
Taq 0.6 $\mu$ l  
W-1 (F1/F2) 0.6 $\mu$ l  
Primer mix 6 $\mu$ l  
H<sub>2</sub>O 93.3  
150 $\mu$ l

For each primer pair use two different DNA samples  
and 1 water control.

## FURTHER OPTIMISATION OF ILIB PCR

$MgCl_2$	1.0 $\mu$ l	1.5 $\mu$ l	5.0 $\mu$ l
Temp.	57 °C	60 °C	63 °C

	<u>MM1</u>	<u>MM2</u>	<u>MM3</u>
$\alpha$ ffor	5 $\mu$ l (45)	5 $\mu$ l (45)	5 $\mu$ l
$MgCl_2$	<del>2.5</del> 1.0 $\mu$ l (9)	1.5 $\mu$ l (13.5)	2.0 $\mu$ l
dNTPs	4 $\mu$ l (36)	4 $\mu$ l (36)	4 $\mu$ l
Tag	0.2 $\mu$ l (1.8)	0.2 $\mu$ l (1.8)	0.2 $\mu$ l
WT	0.2 $\mu$ l (1.8)	0.2 $\mu$ l (1.8)	0.2 $\mu$ l
Primer	5 $\mu$ l (45)	5 $\mu$ l (45)	5 $\mu$ l
$H_2O$	32.6 $\mu$ l (23.4)	32.1 $\mu$ l (28.9)	28.8 $\mu$ l (284.4)
	<u>48 <math>\mu</math>l</u>	<u>48</u>	<u>48 <math>\mu</math>l</u>

### Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: ~~90.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec~~ time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h47m59s

### Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: ~~90.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec~~ time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

### Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: ~~90.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec~~ time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h38m59s

Loading order:

2) A <sub>1</sub> (1)	Temp 57°C	mm 1	DNA sample	96
3) A <sub>1</sub> (2)	"	mm 1	"	97
4) A <sub>1</sub> (H <sub>2</sub> O)	"	mm 1	H <sub>2</sub> O control	
5) A <sub>2</sub> (1)	"	mm 2	"	96
6) A <sub>2</sub> (H <sub>2</sub> O)	"	mm 2	"	97
7) A <sub>2</sub> (H <sub>2</sub> O)	"	mm 3	H <sub>2</sub> O control	
8) A <sub>3</sub> (1)	"	mm 3	"	96
9) A <sub>3</sub> (2)	"	mm 3	"	97
10) A <sub>3</sub> (H <sub>2</sub> O)	"	mm 3	H <sub>2</sub> O control	

H

11) B <sub>1</sub> (1)	Temp 60°	mm 1	# 96
12) B <sub>1</sub> (2)	"	mm 1	# 97
13) B <sub>1</sub> (H <sub>2</sub> O)	"	mm 1	H <sub>2</sub> O control
14) B <sub>2</sub> (1)	"	mm 2	# 96
15) B <sub>2</sub> (2)	"	mm 2	# 97
16) B <sub>2</sub> (H <sub>2</sub> O)	"	mm 2	H <sub>2</sub> O control
17) B <sub>3</sub> (1)	"	mm 3	# 96
18) B <sub>3</sub> (2)	"	mm 3	# 97
19) B <sub>3</sub> (H <sub>2</sub> O)	"	mm 3	H <sub>2</sub> O control


(~~14~~)

20) C <sub>1</sub> (1)	Temp 63°	mm 1	# 96
21) C <sub>1</sub> (2)	"	mm 1	# 97
22) C <sub>1</sub> (H <sub>2</sub> O)	"	mm 1	H <sub>2</sub> O control
23) C <sub>2</sub> (1)	"	mm 2	# 96
24) C <sub>2</sub> (2)	"	mm 2	# 97
25) C <sub>2</sub> (H <sub>2</sub> O)	"	mm 2	H <sub>2</sub> O control
26) C <sub>3</sub> (1)	"	mm 3	# 96
27) C <sub>3</sub> (2)	"	mm 3	# 97
(C <sub>3</sub> (H <sub>2</sub> O))	"	mm 3	H <sub>2</sub> O control

\* ALL PCR REACTIONS ARE BEING CARRIED OUT ON 25μl  
TO MINIMISE WASTAGE OF REAGENTS \*

To ensure that all my reagents are working and I am not doing anything silly I will try out an already proven P.C.R. - (Alisons)

Date	5.9.95.
Number of Samples	3.....
Disease	Any!
P'morph	Taq 1.....
Water	48.8 $\mu$ l
Buffer	8 $\mu$ l
Magnesium	8 $\mu$ l
dNTPs	6.4 $\mu$ l
Primers	8 $\mu$ l
Taq	0.4 $\mu$ l
W-1	0.4 $\mu$ l
Template	2 $\mu$ l



#### 1b-1 B RESULTS Gel loading order P23:-

A number of bands of the correct size were seen  
Those most prominent were:-

- 60°C Annealing - 1.5  $\mu$ l  $MgCl_2^{2+}$

\* - 63°C Annealing - 1.5  $\mu$ l  $MgCl_2^{2+}$

For Gel pictures see above

To ensure that  
I am not doing  
already proven

ing an  
ut an

Date 5.1.91.95.

Buffer ... 3

MgCl<sub>2</sub> ...

dNTP ... Taq ...

W-1 Water 48.8  $\mu$ l

Prim Buffer 8  $\mu$ l

H<sub>2</sub>O Magnesium 8  $\mu$ l


dNTPs 6.4  $\mu$ l

Primers 8  $\mu$ l

Taq 0.4  $\mu$ l

W-1 0.4  $\mu$ l

Template 2  $\mu$ l



### 1.1.1 B RESULTS Gel loading order P23

A number of bands of the correct size were seen. Those most prominent were:

- 60°C Annealing - 1.5  $\mu$ l MgCl<sup>2+</sup>

\* - 63°C Annealing - 1.5  $\mu$ l MgCl<sup>2+</sup>

For Gel pictures see Figure

## OPTIMISATION OF IL12 PCR

Primers F4 / B2

	mm 1	mm 2	mm 3	mm 4	mm 5	mm 6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Mg <sup>2+</sup>	0.5 <sup>4.5</sup>	1.9	1.5 <sup>13.5</sup>	2 <sup>18</sup>	2.5 <sup>22.5</sup>	3 <sup>27</sup>
dNTP's	2	2	2	2	2	2 (18)
Primers	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Taq.	0.1 μl	0.1	0.1	0.1	0.1	0.1 (0.9)
W-1	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
H <sub>2</sub> O	17.3 <sup>155.7</sup>	16.8 <sup>151.2</sup>	16.3 <sup>146.7</sup>	15.8 <sup>142.2</sup>	15.3 <sup>137.7</sup>	14.8 <sup>133.2</sup>
	25	25	25	25	25	25

Predicted annealing Temp: - 53 °C

DNA (Psoriasis) sample: x 6 + H<sub>2</sub>O  
x 6 + H<sub>2</sub>O  
H<sub>2</sub>O

Make 3x amount - for studies at diff Temperatures  
i.e. a total of 9x above

buffer

Carry out above PCR at same Mg<sup>2+</sup> conc  
but at 50 and 56 °C

50 °C set: A<sub>1</sub> = 6 + A H<sub>2</sub>O } labelled 50 °C on side  
B<sub>1</sub> = 6 + B H<sub>2</sub>O }

56 °C set: A<sub>1-6</sub> + A H<sub>2</sub>O } labelled 56 °C on side  
B<sub>1-6</sub> + B H<sub>2</sub>O }

## Summary of 1L1-2 OPTIMISATION

6 different master-mixes @ 6 diff Mg concn:

- 1) 0.5  $\mu$ l
- 2) 1  $\mu$ l
- 3) 1.5  $\mu$ l
- 4) 2.0  $\mu$ l
- 5) 2.5  $\mu$ l
- 6) 3.0  $\mu$ l

⇒ Each master mix was tested out with 2 DNA samples and water controls

⇒ Reactions were carried out at 50-53 and 56 °C

Loading order. IL-1d OPTIMISATION.

1)	Mw Marker	Q		
2)	DNA sample	125	mm 1	50°C
3)	"	135	21	"
4)		125	32	"
5)		135	42	"
6)		125	53	"
7)		135	3	"
8)		125	4	"
9)		135	4	"
10)		125	5	"
11)		135	5	"
12)		125	6	"
13)		135	6	"
14)	H <sub>2</sub> O control	1	3	"
15)	H <sub>2</sub> O control	3	3	"

1)	Mw Marker.			
2)	Sample	125	mm 1	53°C
3)		135	1	"
4)		125	2	"
5)		135	2	"
6)		125	3	"
7)		135	3	"
8)		125	4 - 000	"
9)		135	4	"
10)		125	5	"
11)		135	5	"
12)		125	6	"
13)		135	6	"
14)		125	3	"
15)		135	6	"

1) Molecular wt marker

2) DNA sample 125

Mm 1

56°C

3) " 135

" 1

4) " 125

" 2

5) " 135

" 2

6) " 125

" 3

7) " 135

" 3

8) " 125

" 4

9) " 135

" 4

10) " 125

" 5

11) " 135

" 5

12) " 125

" 6

13) " 135

" 6



14) / " 125

15) / " 135

16) H<sub>2</sub>O control

17) H<sub>2</sub>O control

## RESULTS



at 50°C only smears were observed

Conditions for 1L1D PCR are high temp

low magnesium

53°C 1mM magnesium

1) Molecular wt marker

2) DNA sample 125

mm 1

56°C

3) " 135

" 1

"

4) " 125

" 2

"

5) " 135

" 2

"

6) " 125

" 3

"

7) " 135

" 3

"

8) " 125

" 4

"

9) " 135

" 4

"

10) " 125

11) " 13

12) " 12

13) " 13


14) / " 12

15) / " 13

16) " H<sub>2</sub>O control

17) " H<sub>2</sub>O control

## RESULTS



From Above Gel it seems that best conditions for 1L1 $\times$  PCR are high temp at 50°C only smears were observed

low magnesium

53°C 1mM magnesium

6/9/95

OPTIMISATION OF LIB F2/B1 PRIMERS

Tried already:- 53°C annealing Temp

2.5mM Magnesium

Try now Mg 1 2 3 4 5 6 mM  
Temp 50° 56° 60°

Reagent	1	2	3	4	5	6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5 (22.5)
Mg	0.5	1	1.5	2	2.5	3. (27)
dNTPs	2	2	2	2	2	2. (18)
Primers	2.5	2.5	2.5	2.5	2.5	2.5. (22.5)
Tag	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
W1	0.1	0.1	0.1	0.1	0.1	0.1 (0.9)
H <sub>2</sub> O	16.3 <sup>(146)</sup>	15.8 <sup>(142.2)</sup>	15.3 <sup>(137.4)</sup>	14.8 <sup>(135.2)</sup>	14.3 <sup>(129.7)</sup>	13.8. (124.8)
	24	24	24	24	24	

label samples:-

DNA sample - 157 = ✗  
- 167 = 4

IL-1B OPTIMISATION: Primers F2/B1

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
    step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 2: temp: 50.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 3h 2m23s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
    step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h52m11s

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
    step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 2: temp: 60.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
    step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h46m11s

## Loading Order


	<u>DNA sample</u>	<u>Master mix</u>	<u>Temperature</u>
①	Mw marker		
②	157	mm 1	50°
③	167	1	50°
④	157	2	"
⑤	167	2	"
⑥	157	3	"
⑦	167	3	"
⑧	157	4	"
⑨	167	4	"
⑩	157	5	"
⑪	167	5	"
⑫	157	6	"
⑬	167	6	"
⑭	<del>157</del> H <sub>2</sub> O		"
⑮	<del>167</del> H <sub>2</sub> O		"

①	<del>157</del>		
②	157	mm 1	56°
③	167	1	"
④	157	2	"
⑤	167	2	"
⑥	157	3	"
⑦	167	3	"
⑧	157	4	"
⑨	167	4	"
⑩	157	5	"
⑪	167	5	"
⑫	157	6	"

	<u>DNA Sample</u>	<u>Master Mix</u>	<u>Temperature</u>
⑬	167	6	56°
⑭	+51 H <sub>2</sub> O	6	"
⑮	H <sub>2</sub> O	3	"
⑯			
⑰			
⑱	157 (x)	mm	60°
⑲	167 (x)	1	60°
⑳	157	2	60°
㉑	167	2	60°
㉒	157	3	60°
㉓	167	3	60°
㉔	157	4	60°
㉕	167	4	60°
㉖	157	5	60°
㉗	167	5	60°
㉘	156	6	60°
㉙	167	6	60°
㉚	H <sub>2</sub> O	6	60°
㉛	H <sub>2</sub> O	3	60°
㉜			

- 1<sup>st</sup> gel didn't run properly - will try to Re-run with remaining samples. - 56° low Mg looked good.  
 ⇒ Predicted Annealing Temp = 53°  
 Re-run gel with remainder of samples ⇒ See over.

## RESULTS



IL-1B F2/B1 Primers

OPTIMISED CONDITIONS FOR IL-1B (F2/B1)

= 2 mM magnesium  
56°C annealing Temp.

loading order:-

1)	Q X174 Hae III Mw marker	Temp.
2)	157. (Francis Escorialis) mm1	56°
3)	167. mm1	"
4)	157 " 2	"
5)	167 2	"
6)	157 3	"
7)	167 3	"
8)	157 4	"
9)	167 4	"
10)	157 5	"
11)	167 5	"
12)	157 6	"
13)	167 6	"

IL-1B FURTHER OPTIMISATION. (F. / B.)

⇒ Previously shown that 63°C at 1.0-1.5 mM magnesium produced best results. Will try.

Temperatures: - 62°/64°/66°

MgCl<sub>2</sub> 05/11/1995 1/1.5/2

	MASTER MIX 1	MASTER MIX 2	MASTER MIX 3
BUFFER	2.5 μl	2.5 μl	2.5 μl (15)
Mg	0.15 (3)	0.178 (4.5)	0.181 (6)
dNTP's	2.4 μl (10)	2.4 μl	2.4 μl (12)
primer	2.5 μl	2.5 μl	2.5 μl (15)
Taq	0.1	0.1	0.1 (0.6)
W-1	0.1	0.1	0.1 (0.6)
H <sub>2</sub> O	18.8 (128)	16.05 (96.3)	15.8 (72) 94.8
TOTAL	24 μl	24 μl	24 μl

DNA samples used Psoriasis 164 - P

174 - Q

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 64.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h37m47s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 64.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h40m11s

LOADING ORDER

1)	QX174 Hae III					
2)	Psoriasis sample No. 164	mmi			62°	
3)	" "	174	1		62°	
4)	" "	164	2		62°	
5)	" "	174	2		62°	
6)	" "	164	3		62°	
7)	" "	174	3		62°	
8)	" "	H <sub>2</sub> O control	1		62°	
9)	" "	164	1		64°	
10)	" "	174	1		64°	
11)	" "	164	2		64°	
12)	" "	174	2		64°	
13)	" "	164	3		64°	
14)	" "	174	3		64°	
15)	" "	H <sub>2</sub> O control	2		64°	
16)	" "	164	1		66°	
17)	" "	174	1		66°	
18)	" "	( 164	2		66°	
19)	" "	174	2		66°	
20)	" "	164	3		66°	
21)	" "	174	3		66°	
22)	" "	H <sub>2</sub> O control	3		66°	
23)	QX174 Hae III					

From Yesterdays disastrous cd:-

1) ~~Q~~ x 174 Hae III

2) 157

1 56°

3) 167

1 56°

4) 157

2 56°

5) 167

2 56°

6) 157

3 56°

7) 167

3 56°

8) 157

4 56°

9) 167

4 56°

10) 157

5 56°

11) 167

5 56°

12) 157

6 56°


13) 167

6 56°

14) 157

1 60°


## RESULTS



L1-B F2/B1 Primers

Best band = 2mM Mg

56° Annealing Temp



L1-B F1/B1

Best band = .66 °C

Mg Cl<sup>2+</sup> = 2mM

wrong sized fragment  
must have used

wrong Primers

7/8/95

## IL-12 - F5/B2 OPTIMISATION

Predicted annealing Temperature  $54.6^{\circ}$ 

	Mix 1	Mix 2	Mix 3	Mix 4	Mix 5	Mix 6
Buffer	2.5	2.5	2.5	2.5	2.5	2.5
MgCl <sub>2</sub>	0.5	1	1.5	2	2.5	3.0
dNTP's	2	2	2	2	2	2
Primers	2.5	2.5	2.5	2.5	2.5	2.5
Tag	0.1	0.1	0.1	0.1	0.1	0.1
W-1	0.1	0.1	0.1	0.1	0.1	0.1
H <sub>2</sub> O	16.8	16.05	15.3	14.8	13.8	13.8
	24 $\mu$ l					

Try above mixes at  $52^{\circ}\text{C}$  /  $54^{\circ}\text{C}$  /  $56^{\circ}\text{C}$ .

Make up enough for 8 reactions

	Mix 1	Mix 2	Mix 3	Mix 4	Mix 5	Mix 6
Buffer	20	20	20	20	20	20
MgCl <sub>2</sub>	4	8	12	16	20	24
dNTPs	16	16	16	16	16	16
Primers	20	20	20	20	20	20
Tag	0.8	0.8	0.8	0.8	0.8	0.8
W-1	0.8	0.8	0.8	0.8	0.8	0.8
H <sub>2</sub> O	130.4	128.4	122.4	118.4	114.4	110.4
	192	192	192	192	192	192

## Contents of Step 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 52.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h58m47s

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
~~step 2: temp: 54.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec~~  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h55m47s

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:  
step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
~~step 2: temp: 54.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec~~  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times  
step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h52m11s

loading Order

1)	QX174 marker			
2)	157	1		56°
3)	167	1		56°
4)	157	2		56°
5)	167	2		56°
6)	157	3		56°
7)	167	3		56°
8)	157	4		56°
a)	167	4		56°
10)	157	5		56°
11)	167	5		56°
12)	157	6		56°
13)	167	6		56°
14)				

$\lambda \times 174$  Hae III Molecular Wt marker

1)	157	mm 1	52°
2)	167	mm 1	52°
3)	157	mm 2	52°
4)	167	mm 2	52°
5)	157	mm 3	52°
6)	167	mm 3	52°
7)	167	mm 4	52°
8)	167	mm 4	52°
9)	157	mm 5	52°
10)	167	mm 5	52°
11)	157	mm 6	52°
12)	167	mm 6	52°


H1

$\frac{142}{143}$   
 $\frac{143}{144}$

157	mm 1	54°
167	mm 1	54°
157	mm 2	54°
167	mm 2	54°
157	mm 3	54°
167	mm 3	54°
157	mm 4	54°
167	mm 4	54°
157	mm 5	54°
167	mm 5	54°
157	mm 6	54°
167	mm 6	54°

Results:

54°



From these results

the best bands  
are appearing at low  
Mg Cl<sub>2</sub> and 54°C.

⇒ Try 1mM Mg  
53 & 55°C

+ 2μl DNA sample

← DISASTER  
GEL!!!

IL-1B F1 / B1 OPTIMISATION

Temp  $62^{\circ}$   $64^{\circ}$   $66^{\circ}$

Mg 1 1.5 2


Reagent	mm1	mm2	mm3
Buffer	2.5	2.5	2.5 (20)
Mg	0.5 (4)	0.75 (6)	1 (8)
dNTPs	2	2	2 (16)
Primers	2.5	2.5	2.5 (20)
WT	0.1	0.1	0.1 (0.8)
Tag	0.1	0.1	0.1 (0.8)
H <sub>2</sub> O	18.8 (112.8)	16.05 (96.3)	15.8 (94.8)
	24	24	24

Use 2  $\mu$ l DNA Template.

Make up enough for 8 reactions (quantities in brackets)

loading Order.

1)	$\lambda$ X174 Hae III		62
2)	157	mm 1	62°
3)	167	mm 1	62°
4)	157	mm 2	62°
5)	167	mm 2	62°
6)	157	mm 3	62°
7)	167	mm 3	62°
8)			
8)	157	mm 2	64°
9)	167	mm 2	64°
10)	157	mm 3	66°
11)	167	mm 3	66°
12)	$\lambda$ X174 Mu marker	H1 H2 H3	



1353  
1078  
72

B<sub>1</sub>/F<sub>1</sub>

\* best sample = 66°C Annealing  
2mM Mg

- 803 bp.

It seems that for the larger PCR products using 2 $\mu$ l DNA, works best.

8/9/95 11-12 P.C.R. OPTIMISATION.

F5/B2 (P38 For initial expts).

Try 2mm magnesium at 5B°, 55° and 60°

Reagents	Per reaction	Master mix (12)	
- Buffer	2.5	30	
- MgCl <sub>2</sub>	1.0	24	tubes labelled
- dNTPs	2	24	Gp) R,
- Primers	2.5	30	Ap) S,
- w1	0.1	1.2	
- Taq	0.1	1.2	
- H <sub>2</sub> O	14.8	177.6	
TOTAL.	24.0		

qd 208

F4 | B.2

(INITIAL OPTIMISATION PT 2)

First results not convincing: - Try.

1mM 2mM 3mM MgCl  
Annealing Temp 58° 60° 62°.

REAGENTS	mm1	mm2	mm3
Buffer	2.5	2.5	2.5 (20)
MgCl <sub>2</sub>	0.5 (4)	1.0 (8)	1.5 (12)
dNTPs	2	2	2 (16)
Primers	2.5	2.5	2.5 (20)
W-1	0.1	0.1	0.1 (0.8)
Taq	0.1	0.1	0.1 (0.8)
H <sub>2</sub> O	16.3 (30.4)	15.8 (126.6)	15.3 (122.4)
TOTAL	24.0	24.0	24.0

make enough for 8 reactions

L - 157  
m - 167

14/11/95

2x F1 PCR's using biotinylated primers (for mohammed)

Buff

buffer

21.0  $\mu$ l

MgCl<sub>2</sub>

8.4  $\mu$ l

dNTP

1.68  $\mu$ l

Primer

F1 1.05

W1

0.1 2.1

4.2

Taq

4.2

H<sub>2</sub>O


13.45  $\mu$ l

Do 2 Samples 140 | 188 with parafin oil

- 2 Samples " " with Wax

Running order:-

- 1) QX HacIII
- 2) 140 oil
- 3) 188 oil
- 4) 140 wax
- 5) 188 wax



20/11/95

⇒ PCR up a further (5 x 10g) reactions for two samples. This time to avoid any ethanol ppt steps, the sample will be run straight through a column after band purification - this will miss out any ethanol & also remove excess primer.

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 66.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 2h37m47s

Prepare 2 samples with parafin oil and two with wax, incase parafin oil is interfering with the reaction.

28/11/95

SSCP - Try IL12 - 889 PCR cold to see if it works. - if it does redo hot & follow Tarras protocol

Contents of Step 1:

step 0: temp: 96.0 deg C time: 0h 2m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 45 times:  
step 1: temp: 94.0 deg C time: 0h 16m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 45 times  
step 4: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 53.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 6: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 3h36m 8s

Date...../...../..... 58°C

Number of Samples ..... 10 ~~µl~~

Disease.....

P'morph.....

IL12 - 889 (30µl PCR)

each tube:

3µl 10x Buffer

Water	6.1µl	2.4µl dNTPs
Buffer	1.1µl	1.5µl 603 primer
Magnesium	1.1µl	1.5µl 604 primer
dNTPs.	0.8µl	1.2µl MgCl <sub>2</sub>
Primers M.i.	1.1µl	0.15µl W-1
Taq:	0.05µl	0.15µl <del>K</del> Taq
W-1	0.05µl	18.9µl H <sub>2</sub> O
Template:	.....µl	28.8µl

add 28.8µl to each tube + 1.2µl DNA

30/11/95

Try PCR using Angie / Nicolas Protocol. - and annealing temperature at 50°C

Make up two times the volume - 5 DNA Samples + 1 H<sub>2</sub>O control.

10μl PCR

100

1μl buffer

0.8μl dNTP

0.5μl each Primer <sup>603</sup> <sub>604</sub>

0.4μl Mg<sup>2+</sup>

0.05 μl

0.05 μl Taq

6.3μl H<sub>2</sub>O

Biotinylated primers are correctly biotinylated - START AGAIN

Redo PCR using B primer only

When doing PCR do 45 cycles instead of 35  
this will also help insure against excess primers  
interfering with the reaction.

Do 4 PCR's (5 x 100 μl reactions) 2 normal  
and two cut out of gel & purified

REAGENT	VOLUME	Samples used
Buffer	110 200	
dNTP	88 100	133
primer	55 <del>μl</del> each 100 μl ea.	139
MgCl <sub>2</sub>	44 80	173
Taq	44 8	179
W-1	44 8	
H <sub>2</sub> O	695 <del>μl</del> 1264	
	1100 μl	

— Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 60.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h26m21s

gels PCR

Primer	50 each	- PCR DID NOT
MgCl <sub>2</sub>	40	WORK - TRY
Taq	5	AGAIN WITH
W-1	5	FRESH REAGENTS
H <sub>2</sub> O	650 μl	CSP. dNTP

21/1/96

Re-dissolve <sup>returned</sup> primers in 150µl H<sub>2</sub>O - measure O.D.

5µl → 995µl H<sub>2</sub>O → O.D.

SAMPLE	A320	A280	A260	280/260	260/280	PROTEIN	NUCLEIC ACID
1.0000	0.0036	0.1146	0.1865	0.6065	1.6488	33.653	7.5100
2.0000	0.0048	0.1625	0.2776	0.5782	1.7296	38.189	11.478

B1 - dilution factor =

$$A_{260} \times \text{DILN} \times 0.033/\text{nm} \times 10^6$$

$$= \left( 0.1865 \times 200 \times 0.033 / \frac{7255}{6787.6} \right) 10^6$$

20

1:9 dilution

$$\underline{B2} \quad \left( 0.2776 \times 200 \times 0.033 / \frac{8249}{7275} \right) 10^6$$

20

1:89 dilution

1:11 dilution

Do ~~ILIB~~ PCR's

	m1	m2
Buffer	200	200
MgCl <sub>2</sub>	80	80
dNTPs	160	160
Primer	100 $\mu$ l F2 B1	100 $\mu$ l F5 B2
Taq	8	8
W-1	8	8
H <sub>2</sub> O	1264	1264
	2,000	2,000

4 reactions @ 100  $\mu$ l x 5

Use samples 115 } 1,1 for +3953  
151 } 1,1 for +3953

121 } 2,2 for +3953  
126 }

23/1/96

Running order:-

1) QX Hae III

Run 10  $\mu$ l on

2) 115 F2

1% Agarose gel

3) 151 F2

4) 121 F2

5) 126 F2

6) 115 F5


7) 151 F5

8) 121 F5

9) 126 F5

10) ~~QX~~ H<sub>2</sub>O control

11) QX Hae III



Do 115 and 121 for each Normally

- 151 and 126 extract band from gel to remove primers

### RUNNING ORDER - (ON TAE GEL)

F5 151

F2 151


F2 126

F5 126

- Bands were cut out of the gel, and purified through glass wool and ethanol precipitation

Samples were made single-stranded using dynabeads. Samples of supernatants were kept & run on a 1% Agarose gel.

### RUNNING ORDER



1) QX Hae III

2) 1<sup>st</sup> wash F2 121

3) " " F5 151

4) 1<sup>st</sup> wash SS DNA F2 121

5) SS DNA 121

6) 1<sup>st</sup> wash F2 121

7) " " " 115

8) SS DNA F5 151

We can see from the previous gel picture that for samples which had excess primer removed by running samples down the gel showed no loss of sample on 1<sup>st</sup> wash and a clear band appeared on the lane incorporating ssDNA. - The band may not be very bright but it is well known that ssDNA does not incorporate ssDNA so well.

⇒ Try Sequencing F5 151 DNA.

23/1/96


Sequencing of PCR product (not biotin) is now possible and very successful - prepare PCR's for 10 samples

F<sub>1</sub> - B<sub>1</sub> & F<sub>5</sub> - B<sub>2</sub> 5 Tag 1,1

and 5, X Tag 2,2.

⇒ enough for 150 reactions each

⇒ Reagent	MM1	MM2.
Buffer	500μl	500μl
MgCl <sub>2</sub>	200μl	200μl
dNTP	400μl	400μl
Primer	250μl each	250μl each
W-1	20μl	20μl
Taq	20μl	20μl
H <sub>2</sub> O	2960	2960
	4600	4600



### Running order

- 1) 115
- 2) 121 ✓
- 3) 126 ✓
- 4) 144 ✓
- 5) 165 ✓
- 6) 113
- 7) 139 ✓
- 8) 157 ✓
- 9) 113
- 10) 157 ✓
- 11) H<sub>2</sub>O

#### Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 s

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 s

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 s

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 s


end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 s

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 s

total runtime (approx.): 3h36m 1s

Take all samples that worked, pool their DNA and ethanol extract - then run on a gel and band purify the sample for sequencing



← Run on a gel  
to check samples  
had not been lost  
by above process

1/21/95

MUTS

③

off

96 °C - 5'

37 °C - 5'

X3

- 1) Take 50 $\mu$ l DNA (genomic) [mix 50] x3
- 2) Add 140 $\mu$ l buffer <sup>TRIS HCl pH 7.5</sup> to make volume ~~210~~ 210 $\mu$ l
- 3) Add 10 $\mu$ l beads
- 4) Agitate gently by flicking - incubate @ 37°C - 1/2 hr with shaking
- 5) Centrifuge briefly
- 6) Apply magnet
- 7) Wash in 200 $\mu$ l Wash Buffer X2
- 8) Add buffer from step 2 - heat - 75°C 15' Remove Supernatant - & PCR - Using F1/B1 primers as before! - Run on gel to see if anything's there!



1/2/96

- Do a 100 $\mu$ l PCR reaction with Alison's Taq primers expected fragment size.
- also Do 100 $\mu$ l PCR with ILIB F2-B, expected fragment size - 400 bp

Reagent	Volume	MM
Buffer	12.2	20
Mg <sup>2+</sup>	2	20
dNTP	2.16	16
W <sup>-1</sup>	0.1	1
Tag	0.1	1
Primer	2-mix	20
H <sub>2</sub> O	12.2	122

Reagent	Volume
Buffer	15 $\mu$ l 20ml
Mg <sup>2+</sup>	6ml 8ml
dNTP	16 $\mu$ l
Primer	10 $\mu$ l each
Tag	0.8 ml
W <sup>-1</sup>	0.8 ml
H <sub>2</sub> O	120 $\mu$ l

10ul template

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h36m 1s

5/2/96

Samples for sequencing:- 1 - 121

2 - 126

3 - 157

4 - 151

5 - 144

6 - 139

7 - 165

- Give Hazel F2 primer to sequence half of each of the samples: - Tag FS sequencing

Do PCR's as before ←

Contents of Sto 2:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 94.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h36m 1s

Contents of Sto 1:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

begin cycle 45 times:

step 1: temp: 96.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 2: temp: 64.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec

end cycle 45 times

step 4: temp: 72.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec

step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec

total runtime (approx.): 3h24m27s

Do PCR's on  $\beta$  for direct sequencing from PCR product.

Do PCR with  $F5/B2$  and sequencing using both  $F6/F4$  primers

$F1/B1$  PCR - again both  $F1/F2$  primers

REAGENT	Volume	Volume
10 reactions $\times 10$		
buffer	500 $\mu$ l	500
= 1000	200 $\mu$ l	400
$MgCl_2$	200	200
Primer	250 each	250 each
Tag	20 $\mu$ l	20 $\mu$ l
W-1	20 $\mu$ l	20 $\mu$ l
$H_2O$	2960 $\mu$ l	2960 $\mu$ l
	4600 $\mu$ l	4600 $\mu$ l

Running order:-

1) QX Hoe III

2) 115 F

3) 151

4) 157

5) 165

6) 174

7) 121


8) 126

9) 189

10) 144

11) 152

12)



3/2/96

- Run PCR products on gel - band extract and  
ethanol precipitate - Give to Hazel to  
sequence

Gel to check products are still there

### Running Order

1) QX Hae III marker

2) F5 174

3) F1 126

4) F5 151

5) F1 151

6) F1 176

7) F5 152

8) F5 157

9) F1 139

10) F5 121

11) F5 139

12) F5 165

13) F1 152

14) F1 165


15) F5 115

16) F5 126

17) F1 157

18) F1 144

19) F1 115




10/2/96

Do Biotinylated PCR's for 5x 1,1  
5x 2,2.

<u>REAGENT</u>	<u>VOLUME</u>
buffer	500 $\mu$ l
Mg Cl <sub>2</sub>	200 $\mu$ l
dNTP	400 $\mu$ l
Primer	250 $\mu$ l each (used 200 F1 100 BBI)
W-1	20 $\mu$ l
Taq	20 $\mu$ l
H <sub>2</sub> O	2960 $\mu$ l

Run samples on a 1% agarose gel



RUNNING ORDER

- 1) QX Hinf
- 2) 115 F1 | B61
- 3) 151 F1 | B61
- 4) 157 F1 | B61
- 5) 165 F1 | B61
- 6) 174 F1 | B61
- 7) 121 F1 | B61
- 8) 126 F1 | B61
- 9) 139 F1 | B61
- 10) 144 F1 | B61
- 11) 152 F1 | B61

11/12/96

Pool all 5 x 100  $\mu$ l's and ethanol precipitate resuspend in 20  $\mu$ l for running ~~of~~ on TAE gel & band extract

Running Order

- 1) 139
- 2) 151
- 3) 126
- 4) 144
- 5) 115
- 6) 121
- 7) 157
- 8) 174
- 9) 152
- 10) 165

After band extraction, samples were ethanol precipitated and resuspended in 40  $\mu$ l H<sub>2</sub>O. 1  $\mu$ l of sample was run on a 1% agarose gel to check that I have not lost them during ethanol precipitation.



4/2/96

Do 1L-1 Tag PCR on 50 patients ( $4 \times 100\mu\text{l}$ )  
for control experiments using mab 5 beads  
⇒ Worked OK.

Mohammeds Solutions -F1 LS7  
-F1-174  
F1-152  
F1-121  
-F1-139

$$1\% = 1\text{ml} / 100$$

$$0.005\% =$$

$$0.005\%$$

$$\text{mab } 5\mu\text{l} \rightarrow 100\text{mls}$$

$$0.5 \rightarrow 10\text{mls}$$

$$20\mu\text{l} \rightarrow 10\text{mls}$$

$$1\mu\text{l} \rightarrow 40$$

—

15/2/96

## MUTS - New Protocol

- Carry out reaction in 10mM Tris-HCl (Not 1M)
- Put  $^{0.005\%}$  Tween 20 in buffer rather than BSA
  - ensure buffer is FILTERED through 0.45μM
- Binding occurs at room temp - Not 37°C
- Wash beads with 10mM Tris-HCl (No Tween)

### CONTROLS

(1) Beads - No DNA - boil, remove 5' - PCR as a H<sub>2</sub>O control.

(2) One @ 37°C and one at Room temp

beads with no BSA but Tween were OK

### REACTION BUFFER

- 1M Tris HCl pH 7.5 }  
- 250mM MgCl<sub>2</sub> } 50X       $\hookrightarrow$  PCR Buffer!  
- 5mM dTT }  
- 0.5mM EDTA }

+ add MgCl<sub>2</sub>

~~After~~ Equal Volume

Block beads with 5mM P.V.P.

2g  $\rightarrow$  10mls


~~samples~~  
- Again samples were heated to 75°C for fifteen minutes and supernatant removed immediately

2/1/2/96


Sequencing has revealed a possible G-C change.

Design method for screening possible polymorphs


- No enzyme (except <sup>cont by</sup> Tsp RI) cuts at the altered site. Therefore a cut site must be engineered - one possibility is to engineer a cut site for *Hinf* I in a primer. There are 2 intrinsic *Hinf* I sites + the engineered one.



#### EXPECTED BANDING PATTERNS



Design oligo which creates the *Hinf* site.  
Because,



21/2/96

Column #3 MM 828

Seq: Seq03

Overall: 50.5

ASWY: 98.0

Num Base ASWY

2 C	100.0
3 C	98.5
4 A	96.0
5 T	97.0
6 T	97.6
7 T	98.0
8 A	98.0
9 A	98.0
10 A	97.9
11 T	98.0
12 A	98.1
13 A	97.4
14 T	97.6
15 C	97.8
16 T	97.9
17 G	98.0
18 A	98.1
19 G	98.2
20 C	98.3
21 T	98.2
22 T	97.4
23 A	97.5
24 T	97.7
25 A	97.8
26 T	97.8
27 C	97.9
28 G	98.0
29 T	98.1
30 A	98.1
31 T	97.7
32 G	97.8
33 A	97.9
34 G	98.0
35 T	98.0

SEQUENCE: t6911 Primer

5' TGA GTT TTA TAT ATT ATT  
CGA GTC TAA TAA ATT TAC C 3'

Predicted annealing temperature  
57°C

(Using hypercard)

Column 1

14:53:21 , 23/ 2/96

Run ID : MW 831  
Cycle : 002 UMOL  
End Proc: End CESS (DMT = Off)  
Sequence: Seq01

Carol

Average

Step-wise

Yield : 98.8

Total bases = 35

A= 11, G= 4, C= 5, T= 15, 5= 0, 6= 0, 7= 0, 8= 0  
(mixed bases= 0)

MW: 10706.0

5' > CCC ATT TAA ATA ATC TGA GCT TAT ATA TTT TGA GT <3'

Ethanol precipitate primers and set  
up optimisation of PCR's.

#### OPTIMISATION STRATEGY

	2mM	3mM	4mM	5mM
Buffer	35	35	35	35
MgCl <sub>2</sub>	14 <sub>u</sub> l	21 <sub>u</sub> l	28 <sub>u</sub> l	35 <sub>u</sub> l
dNTP	28 <sub>u</sub> l	28 <sub>u</sub> l	28 <sub>u</sub> l	28 <sub>u</sub> l
Primer	35 <sub>u</sub> l	35 <sub>u</sub> l	35 <sub>u</sub> l	35 <sub>u</sub> l
W-1	1·4	1·4	1·4	1·4
Taa	1·4	1·4	1·4	1·4
H <sub>2</sub> O	207·2	207·2	207·2	207·2

2 Samples - 3 temperatures + 1 h<sub>2</sub>O  
= 7 tubes

$$A_{260} = 0.2411$$

$$(A_{260} \times \text{dil}^2 \times 0.033 / \text{MW} \times 10^3) = \text{dilution factor}$$

$$\left( \frac{0.2411 \times 200 \times 0.033}{107.06} \right) \times 10^3 = \text{Dilution factor}$$

1.59 mg/ml


0.148 mm

1.1.7.43


### RUNNING ORDER:

- 1) QX Hae III
- 2) 255 } 1 mM MgCl<sub>2</sub> 54°C
- 3) 256 } 4 mM MgCl<sub>2</sub>
- 4) 255 } 1 mM Mg 56°C
- 5) 256 } 4 mM MgCl<sub>2</sub>
- 6) 255 } 1 → 4 mM MgCl<sub>2</sub>
- 7) 250 } - 58°C
- 8) H<sub>2</sub>O - 1
- 9) - 2
- 10) - 3
- 11) - 4

Run 10<sub>ul</sub> on a  
1% Agarose gel



- See p129 for  
Running order



all Magnesium concentrations noted at all  
 3 temperatures Do  $30 \times 50\mu\text{l}$  PCR's  
 on  $10 \times 1,1 \}$   
 $10 \times 1,2 \}$  for Taq Polymorphism  
 $10 \times 2,2 \}$   
 i.e. 31 tubes

Reagent	Volume	
Buffer	155 $\mu\text{l}$	
MgCl <sub>2</sub>	62 $\mu\text{l}$	
dNTP	124 $\mu\text{l}$	1 $\mu\text{l}$ template
Primer	77.5 each	+ 4.6 $\mu\text{l}$ master mix
W-1	6.2 $\mu\text{l}$	
Taq	6.2 $\mu\text{l}$	Carry out reaction
H <sub>2</sub> O	917.6	@ 56 °C
	1426	

Digestions  $2.5 \mu\text{l}$  enzyme per tube  
 $-3060$  tubes

7.5  $\mu\text{l}$  enzyme  $2.25 \mu\text{l}$  15  $\mu\text{l}$  enzyme  
 $60 \mu\text{l}$  buffer }  $120 \mu\text{l}$  buffer }  
 20 ml

leave to digest overnight @ 37 °C

Run samples on a 9% polyacrylamide gel  
 @ 200V for 30 minutes

Samples used: -

TAQ + 3845

Gel 1

112 - 1,1 - CC

113 - 1,1 - -

114 - 2,2 - GG

115 - 1,1 - CC

116 - 1,2 - GC

117 - 1,1 - CC

118 - 1,2 - CG

119 - 1,1 - -

120 - 1,1 - CC

121 - 2,2 - GG

122 - 1,2 - GC

123 - 1,2 - GC

124 - 1,2 - GC

125 - 1,2 - GC

126 - 2,2 - GG

Gel 2

127 - 1,1 - CC

128 - 1,2 - GC

129 - 1,1 - CC

130 - 1,2 - GC

132 - - CC

133 - - GC

134 - 1,2 - GG

135 - 1,2 - GC

136 - 1,1 - CC

137 - 1,2 - GC

138 - 1,1 - CC

139 - 2,2 - GG

140 - 1,1 - -

141 - 1,1 - -


142 - - -

12 x 1,1


11 x 1,2

4 x 2,2

Thus it appears that Allele (G) of my polymorphism is 100% associated with allele 2 of alisons Taq - and C is 100% linked to allele 1.



-35



## FUTURE PLANS FOR +3845 Polymorphism

- Fusions of AU Rich region to B globin

Normal mRNA  $1/2$  life  
= 17 hrs

- Band shift assays with both mRNAs  
incubated with nuclear extract

1<sup>st</sup> - ASTA using Adcel's Taq-man System


1/3/96

Set up Tag PCR again

6x2x100  $\mu$ l PCR's  $\rightarrow$  4 mixtures  
Reactions

REAGENT	VOLUME
Buffer	40
MgCl <sub>2</sub>	40
dNTP	32
w-1	2
Tag	2
Primer	40
H <sub>2</sub> O	244

- add 8  $\mu$ l template  
for each:-



- Obviously, this has  
not worked - Why am  
I getting such smears

#### CONTROL EXPTS

- 1) Band extract DNA before putting over the column
- 2) Try Re-PCR'ing the PCR product to see if I get the smears then
- 3) Run on acrylamide & stain for protein

8845  
1941  
8904  
20 each

Plan Expts


- PCR - up 1,1 )

2,2 } check that  
1,2 } they are same for  
more

- How much blood?
- Time Scale?

1C1B peaks @ 6 hr

+ 18-24 hr



Tagman 5/10/15 x 20 cycles - Standard curve  
Quantitation of steady state mRNA levels  
⇒ Stability or accumulation of mRNA.

eosinophils - acidic - stain eosin 165

basophils - basic


neutrophils

monocyte - nucleus almost complete cell division  
lymphocyte - from lymph - dense chromatin  
monocyte - kidney-shaped nucleus  
- much less dense than lymphocyte  
- because DNA is ~~not~~ dense

100 cells - 30 neutro  
50 Lympho  
mono

5/3/96


Do Tag PCR for further M13-5 expts  
- as previously set out (p115)



RUNNING ORDER

- 1) QX 1174 Hae III
- 2) Tag
- 3) Tag
- 4) H<sub>2</sub>O control

Also F5-3G<sub>2</sub> PCR for sequencing of 1612 3'



Running Order

~~QX Hae III~~

- 2) 165
- 3) 174
- 4) 149
- 5) 115
- 6) 157
- 7) 152
- 8) 139
- 9) 126
- 10) 121
- 11) 144
12. H<sub>2</sub>O

DO 40 x 25  $\mu$ l reactions to confirm linkage of + ?? to Tag polymorphism

### Tag PCR

3 samples from muts

1 positive controls  $\Rightarrow$  20  $\mu$ l reaction

1 H<sub>2</sub>O control

REAGENT	VOLUME
Buffer	8
MgCl <sub>2</sub>	8
dNTP	12.8
Primer	8 $\mu$ l
Tag	0.8
W-1	0.8
H <sub>2</sub> O	97.8

### RUNNING

- 1)  $\alpha$ x Hae III
- 2) PVP 50
- 3) PVP 100
- 4) Tween
- 5) +ve control
- 6) H<sub>2</sub>O control

Do <sup>48</sup> ~~96~~ x 25  $\mu$ l PCRs of new Polymorphism PCR

REAGENT	volume
Buffer	120
MgCl <sub>2</sub>	96 (48)
dNTP	96
Primer	4120
Taq	48
WT	48
H <sub>2</sub> O	710.4
	110.4

2  $\mu$ l template per reaction

\* Using the plate, some of the samples have evaporated!

- For digestion do 1 $\mu$ l ~~WT~~ 10  $\mu$ l of sample digested with 0.5  $\mu$ l enzyme overnight and a further 1 $\mu$ l with 0.25  $\mu$ l enzyme overnight

mm1

6  $\mu$ l enzyme

18  $\mu$ l buffer

mm2

3  $\mu$ l enzyme

18  $\mu$ l buffer

7-19 0.5  $\mu$ l enzyme / 10  $\mu$ l

21-36 0.25  $\mu$ l enzyme / 10  $\mu$ l

6/3/96

Re-do F5-B62 PCR's on samples:-

152

139

126

121

144

REAGENT	VOLUME
Buffer	130
MgCl <sub>2</sub>	52
dNTP	104
Primer	6.5 $\mu$ l each 6.5 $\mu$ l each
W-1	5.3 $\mu$ l
Taq	5.2 $\mu$ l
H <sub>2</sub> O	769.6
TOTAL	1196

RUNNING ORDER

1) QX HaeIII

2) 126

3) 139


4) 144

5) 152

6) 121

7) H<sub>2</sub>O

8) 115 SS



7/3/96

Do 2 10 + 1H<sub>2</sub>O : 50μl PCRS for new polymorphism:-


REAGENT	VOLUME
Buffer	55 μl
MgCl <sub>2</sub>	22 μl
dNTP	44 μl
Primer W-1	55 μl 27 μl each
Taq	2.2
H <sub>2</sub> O	330 μl 4 μl DNA per react

Contents of Sto 3:

step 0: temp: 96.0 deg C time: 0h 5m 0s slope: 2.00 deg/sec time-inc: 0 sec  
begin cycle 35 times:

step 1: temp: 95.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 2: temp: 56.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 3: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
end cycle 35 times

step 4: temp: 72.0 deg C time: 0h 1m 0s slope: 2.00 deg/sec time-inc: 0 sec  
step 5: temp: 4.0 deg C time: infinite slope: 2.00 deg/sec time-inc: 0 sec  
total runtime (approx.): 2h49m21s



Set up digests o/n.


5 - 0.5  $\mu$ l enzyme / 10  $\mu$ l


5 - 0.25  $\mu$ l enzyme / 10  $\mu$ l

Reaction mix 1: 5  $\mu$ l enzyme + 15  $\mu$ l buffer

2: 2.5  $\mu$ l enzyme + 15  $\mu$ l buffer

leave to cleave @ 37°C overnight  $\Rightarrow$  Run on  
9% polyacrylamide gel.





Unable to cleave the top band even with 10ul PCR product and 0.5μl enzyme @ 37°C overnight.

all samples run today were in equilibrium with +3953 polymorphisms!

7 - 22 - GG

8 - 11 - CC

9 - 1,2 - GC

10 - 1,1 - CC

12 - 1,1 - CC

13 - 1,2 - GC

14 - 1,1 - CC


15 - 1,2 - GC

16 - 1,2 - GC

13/3

OPTIMISATION OF SCREENING FOR +6911

- A dilute 7μl PCR product to 30μl with H<sub>2</sub>O to cut out effects of buffers in PCR on enzyme.
- Add 3.5μl of buffer to each tube
- + 2μl enzyme. leave @ 37°C o/n.



Run 1/2 leave other 1/2 overnight for another night!



134

18/3/95



Do PCR's on 10x F. B.

Buf Reagent	Volume
BUFFER	500
dNTP	400
Mg <sup>2+</sup>	200
Primer	250 each - 100 $\mu$ l biotinylated
Taq	20
W1	20
H <sub>2</sub> O	2960
	4600 $\mu$ l


Redo Tag

P3/82


1) Running order  
QX haem  
2) 153  
3) 126  
4) 114  
5) 144  
6) 139  
7) 121



SQUIGGLES of: 111bb1t connect March 18, 1998 17:07 E9-U  
FOLDRNA of: 111bb1t Check: 877 from: 1 to: 871 March 7, 1998 17:34  
Length: 871 Energy: -229.5





6



## MUTANT

136.

SOURCELES of: Illibbit.connect March 18, 1998 20:57:5-5  
FOLDRNA of: Illibbit Check: 708 from: 1 to: 871 March 18, 1998 20:12  
Length: 871 Energy: -228.8



WT

101

22/3/96

Take Martin's IL-1B probe in PUC 18 plasmid  
and plate out on Amp/met + IPTG + X-gal. Grow  
@ 37°C overnight - Take one white colony  
& inoculate 3 ml LB culture

Ampicillin

Meth

### LB Recipe 1 litre

10g Bacto tryptone

5g Yeast extract

5g NaCl

Autoclave before use.

3 ml culture was used to inoculate  
a large 1 litre culture (2 x 500 ml). Again  
ampicillin and methicillin were used.

Plasmid DNA is extracted by PEG  
precipitation to give good, pure yields.

Using Ruth Herbst's Protocol.

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- BLACK BORDERS**
- IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- FADED TEXT OR DRAWING**
- BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- SKEWED/SLANTED IMAGES**
- COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- GRAY SCALE DOCUMENTS**
- LINES OR MARKS ON ORIGINAL DOCUMENT**
- REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- OTHER: \_\_\_\_\_**

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**